

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 007093	Instrument Location					
Reason for Install/Repair						
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change						
A Other (Specify) annual Inspection at Crime Laboratory						
Check When Done:	U					
1. Surge Protector Installed/Property Grounded.						
✓ 2. Telephone Line Connected to Intoxilyzer® 8000. Fax						
☑ 3. Breath Tube Heated.						
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).						
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).						
☑ 6. Run Tests:						
A. Print Test (Level1, Function P).						
B. ACA Test (Level 1, Function C).						
A N K. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).				
7. Repair and/or Maintenance Performed (if any):	val Inspection					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5		ntoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	G) at the Intoxilyzer® Location at the Age	ncy.				
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form				
Field Inspector Signature		Date 4/29/2024				
Crime Laboratory Use Only						
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on for County of Burleigh, North Dakota, is certified to be a true and correct co	ile at the Office of Attorney General, Crin					
Reviewed/Certified By		Certified Date 30April2024				

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007093
Location = BOPD 8164.16.00 09/18
04/29/2024 17:55

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 14323080A4 Standard Cyl #? Standard Expiration? 06/05/2025 Oper No? 130538

Flow Cal. Date: 06/07/2019 Slope 638 Intercept -611758

IR Calibration Date: 12/16/2019

3um 9um

Oth Coef(*100): -14610 -32356

1st Coef(*100): 264615 138933

2nd Coef(*100): 2992 1364

H2O adj(mq/l*10k): 369 621

****** Printer Test End ********

Operator Signature
TREVOR BERGERSON

Remarks:

Install after annual inspection.

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BOPD
04/29/2024

SN 80-007093 8164.16.00 09/18 17:55

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	17:55
02	Std.	Gas	0.081	17:56
03	Room	Air	0.000	17:56
04	Std.	Gas	0.082	17:57
05	Room	Air	0.000	17:57
06	Std.	Gas	0.081	17:58
07	Room	Air	0.000	17:58

Lot No = 14323080A4

Cyl No = 2

Exp Date = 06/05/2025

County = 06

Oper No. = 130538

Operator Signature
TREVOR BERGERSON

Remarks:

Install after annual inspection.

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BOPD
04/29/2024

Alcohol Analyzer SN 80-007093 8164.16.00 09/18 17:59

Test			AC	Time
01	Room	Air	RFI*	18:01
02	Room	Air	0.000	18:01

*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI TEST

Sub DOB = 06/03/1988

Sub Sex = Unknown

Weight = NA

Test = OTH

Cit = NA

Dr. Lic. = ND/NA

Lot No = 14323080A4

Cyl No = 2

Expiration Date = 06/05/2025

County = 06

Oper No. = 130538

Operator Signature TREVOR BERGERSON

Remarks:

Install after annual inspection.

Form 106-I8000