



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Kari Stuart</i>		
Location <i>McKenzie County Correctional Facility-WCPD</i>	Intoxilyzer® Serial Number <i>20-004942</i>	
Gas Lot Number <i>14323080AA</i>	Gas Cylinder Number <i>12</i>	Gas Expiration Date <i>00/05/25</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. *0.078* AC  
 B. *0.078* AC  
 C. *0.070* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>04/03/24</i>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <i>08May2024</i>

*-did not have correct location code. Reached out and had correct location code added and a new print test sent. Print test is attached. -AEN*


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004948  
Location = TOXL      8164.14.00 09/16  
04/03/2024      06:05

DRY CAL CHECK


Test	AC	Time
01 Room Air	0.000	06:05
02 Std. Gas	0.078	06:06
03 Room Air	0.000	06:06
04 Std. Gas	0.078	06:07
05 Room Air	0.000	06:07
06 Std. Gas	0.078	06:07
07 Room Air	0.000	06:08

Lot No = 14323080A4  
Cyl No = 18  
Exp Date = 06/05/2025  
County = 27      Oper No. = 132007

  
Operator Signature  
KARI STUART

Remarks: ACA

Form 106-I8000

STATE OF NORTH DAKOTA } SS  
COUNTY OF MCKENZIE  
I certify that the foregoing is a true and correct  
copy of the original document on file in this office.  
Dated at Watford City, North Dakota.  
on this 3rd day of April 2024  
  
Watford City Police Department

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004948  
Location = WCPD                    8164.14.00 09/16  
04/09/2024                            18:26

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            14323080A4  
Standard Cyl #?            18  
Standard Expiration?       06/05/2025  
Oper No?                    133050

Flow Cal. Date:             08/28/2018  
Slope                        642  
Intercept                    -606282

IR Calibration Date:        04/01/2024  
                                    3um                            9um

	3um	9um
0th Coef(*100):	-23587	-22037
1st Coef(*100):	256489	131821
2nd Coef(*100):	3363	1546
H2O adj(mg/l*10k):	339	298

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
MATTHEW EARL

Remarks: *Print test*