

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	
80-006673 Reason for Install/Repair	Both head to Jail	
Reason for install/Repair		
☐ Install After Receiving From Crime Laboratory ☐ Install After		
Other (Specify) Install after Annual In	spection	
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).	
☑ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
A. Print Test (Level1, Function P).		
☑ B. ACA Test (Level 1, Function C).		
	Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):	A	
⊠ 8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.
☑ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G)	6) at the Intoxilyzer® Location at the Age	ncy.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature Kento	n Kussu-	Date 9/30/2024
Crime Laboratory Use Only		,
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on ficulty of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorney General, Crin	
Reviewed/Certified By		Certified Date 01May2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006673
Location = BOTT 8164.14.00 09/16
04/30/2024 14:54

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled DABACABA Start Test Sequence: Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No No Display Volume? Disable On Memfull? Yes 1 # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 01923080A3 Standard Lot #? Standard Cyl #? 052 Standard Expiration? 02/05/2025 010907 Oper No?

Flow Cal. Date: 03/30/2017 Slope 697 Intercept -667725

IR Calibration Date	04/15/20 3um	024 9um
Oth Coef(*100): 1st Coef(*100): 2nd Coef(*100): H20 adj(mg/l*10k):	-14904 290754 2231 289	-25220 143378 1156 324

******* Printer Test End ********

Operator Signature KENTON KOSSAN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006673 Location = BOTT 8164.14.00 09/16 04/30/2024

15:04

DRY CAL CHECK

Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.079 0.000 0.079 0.000 0.078 0.000	15:04 15:05 15:05 15:06 15:06 15:07
0 / 2.0 0 // 2.00		

Lot No = 01923080A3

Cyl No = 052

Exp Date = 02/05/2025

County = 05

Oper No. = 010907

Operator Signature KENTON KOSSAN

Keen m

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006673
Location = BOTT 8164.14.00 09/16
04/30/2024 15:08

Test	AC	Time
01 Room Air	RFI*	15:09
02 Room Air	RFI*	15:09

*Invalid Test Inhibited - RFI

Sub Name = CHRIST, JESUS H Sub DOB = 01/01/1900

Sub DOB = 01/01/1900 Sub Sex = Male Test = OTH

Weight = 150 Cit = 00000001

Dr. Lic. = ND/CHR000001 Lot No = 01923080A3

Cyl No = 052

Expiration Date = 02/05/2025

County = 05

Oper No. = 010907

Operator Signature KENTON KOSSAN

Kern Ku

Remarks:

Form 106-I8000