



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Thomas, Patrick</i>		
Location <i>Fargo Police Department</i>	Intoxilyzer® Serial Number <i>80-004942</i>	
Gas Lot Number <i>01923080A3</i>	Gas Cylinder Number <i>23</i>	Gas Expiration Date <i>2/5/2024</i> <i>2025 aen</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.081 AC
 - B. 0.080 AC
 - C. 0.080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Patrick Thomas</i>	Date <i>5/1/2024</i>
Reviewed By (Crime Laboratory Use Only) <i>A. En</i>	Date <i>02May24</i>

- incorrect location code. Informed officer and requested a print test with the corrected location code. Print test attached. -AEN 02May24


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = FPD 8164.14.00 09/16
05/01/2024 02:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	02:47
02 Std. Gas	0.081	02:47
03 Room Air	0.000	02:48
04 Std. Gas	0.080	02:48
05 Room Air	0.000	02:49
06 Std. Gas	0.080	02:49
07 Room Air	0.000	02:49

Lot No = 01923080A3
Cyl No = 23
Exp Date = 02/05/2025
County = 09 Oper No. = 130056

 #499

Operator Signature
PATRICK THOMAS

Remarks: *Good test, within range.*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = FAPD 8164.14.00 09/16
05/01/2024 09:36

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 01923080A3
Standard Cyl #? 23
Standard Expiration? 02/05/2025
Oper No? 130056

Flow Cal. Date: 05/25/2018
Slope 626
Intercept -548004

IR Calibration Date: 04/09/2024
 3um 9um

0th Coef(*100): -22609 -23575
1st Coef(*100): 260445 132718
2nd Coef(*100): 2885 1338
H2O adj(mg/l*10k): 492 488

***** Printer Test End *****

 #499

Operator Signature
PATRICK THOMAS

Remarks:

Form 106-I8000