

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

	1			
Serial Number 80 -004942	Instrument Location			
Reason for Install/Repair	Fargo Pelice Department - 105 25 th St N, FGO			
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change			
Other (Specify)	or account on ange			
Check When Done:				
☑ 1. Surge Protector Installed/Property Grounded.				
☑ 2. Telephone Line Connected to Intoxilyzer® 8000.				
☑ 3. Breath Tube Heated.				
☑ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio During Test).			
	tenance performed. Install after move from Cass Jail.			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G				
10. Send the Following to the Crime Laboratory: Completed Into	oxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form			
104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature	Date			
Fith The	5/1/2024			
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be	sused for the analysis of broath to determine cleabel			
the date the held inspector performed the installation. This record on file	e at the Office of Attorney General Crimo Laborator, Division in the			
County of Burleigh, North Dakota, is certified to be a true and correct cop	by of the documents received.			
Reviewed/Certified By	Certified Date			
Com Mill	OaMaya4			
	V			

-incorrect location code. Informed officer and requested a print test with the correct location code. Print test attached.

FORM 104-G

-AEN DAMay24

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = FPD 8164.14.00 09/16
05/01/2024 02:41

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Dry Select Std (D/W/I)? Standard Value? 0.080 Standard Lot #? 01923080A3 Standard Cyl #? 23 Standard Expiration? 02/05/2025 Oper No? 130056

Flow Cal. Date: 05/25/2018 Slope 626 Intercept -548004

IR Calibration Date: 04/09/2024
3um 9um

Oth Coef(*100): -22609 -23575
1st Coef(*100): 260445 132718
2nd Coef(*100): 2885 1338
H2O adj(mg/l*10k): 492 488

****** Printer Test End *******

Operator Signature PATRICK THOMAS

Remarks: Good test.

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = FPD
05/01/2024

Alcohol Analyzer SN 80-004942 8164.14.00 09/16 02:46

DRY CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	02:47
02	Std.	Gas	0.081	02:47
03	Room	Air	0.000	02:48
04	Std.	Gas	0.080	02:48
05	Room	Air	0.000	02:49
06	Std.	Gas	0.080	02:49
07	Room	Air	0.000	02:49

Lot No = 01923080A3

Cyl No = 23

Exp Date = 02/05/2025

County = 09

Oper No. = 130056

Operator Signature PATRICK THOMAS

Remarks: Good test, within range.

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = FPD
05/01/2024

Alcohol Analyzer SN 80-004942 8164.14.00 09/16 03:38

Test	AC	Time
01 Room	RFI*	03:40
02 Room	0.000	03:41

*Invalid Test Inhibited - RFI

Sub Name = INTSTALLATION, TEST NA

Sub DOB = 01/01/1901

Sub Sex = Male

Weight = NA

Cit = NA

Test = DUI Dr. Lic. = ND/NA

Lot No = 01923080A3

Cyl No = 23

Expiration Date = 02/05/2025

County = 09

Oper No. = 130056

Operator Signature PATRICK THOMAS

Remarks: RFI test Successful

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = FAPD 8164.14.00 09/16
05/01/2024 09:36

******* Printer Test *******

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Display Prelim No. 1.
Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01923080A3 Standard Cyl #? 23 Standard Expiration? 02/05/2025 Oper No? 130056

Flow Cal. Date: 05/25/2018 Slope 626 Intercept -548004

IR Calibration Date: 04/09/2024 3um 9um

Oth	Coef(*100):	-22609	-23575
1st	Coef(*100):	260445	132718
2nd	Coef(*100):	2885	1338
H20	adj($mg/l*10k$):	492	488

******* Printer Test End *******

Operator Signature
PATRICK THOMAS

Remarks: