

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-007098	Instrument Location Tipas PD (Booking)
Reason for Install/Repair	Ting PD (Booking)
│ │ │ │ │ │ │ │ │ │ │ │ │ │ │ │ │ │ │	fter Location Change
	- Crime Laboratory
Check When Done:	O
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Local	ation; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode o	or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any): Non-e	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-6	G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	ntoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
Field Inspector Signature	Date
John helse	5/9/2024
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct county.	
Reviewed/Certified By	Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007098
Location = TIOG 8164.16.00 09/18
05/09/2024 15:46

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() \_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Yes Display Prelim Rslt? Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? NoDisable On Memfull? Yes # of Print Copies? Dry Select Std (D/W/I)? 0.080 Standard Value? Standard Lot #? 28423080A3 Standard Cyl #? 46 Standard Expiration? 11/05/2025 130977 Oper No?

Flow Cal. Date: 09/26/2018 Slope 684 Intercept -588446

IR Calibration Date: 09/26/2018
3um 9um

Oth Coef(\*100): -40948 -15551
1st Coef(\*100): 237555 134140
2nd Coef(\*100): 2571 1433
H2O adj(mg/l\*10k): 1097 457

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature JOSHUA NELSON

Remarks:

I certify this is a true and accurate copy.

Form 106-I8000

X Lande

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TIOG
05/09/2024

SN 80-007098 8164.16.00 09/18 15:51

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:51
02 Std. Gas	0.081	15:52
03 Room Air	0.000	15:52
04 Std. Gas	0.080	15:53
05 Room Air	0.000	15:53
06 Std. Gas	0.081	15:53
07 Room Air	0.000	15:54

Lot No = 28423080A3

Cyl No = 46

Exp Date = 11/05/2025

County = 53

Oper No. = 130977

Operator Signature JOSHUA NELSON

Remarks:

I certify this is a true and accurate copy.

Form 106-I8000

× free

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TIOG
05/09/2024

Alcohol Analyzer SN 80-007098 8164.16.00 09/18 15:47

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:48
02 Std. Gas	0.081	15:48
03 Room Air	0.000	15:48
04 Std. Gas	0.081	15:49
05 Room Air	0.000	15:49
06 Std. Gas	0.081	15:50
07 Room Air	RFI*	15:50
08 Room Air	0.000	15:50

\*Invalid Test Inhibited - RFI

Lot No = 28423080A3

Cyl No = 46

Exp Date = 11/05/2025

County = 53

Oper No. = 130977

Operator Signature JOSHUA NELSON

Remarks:

I certify this is a true and accurate copy.

Form 106-I8000

RFI fifrom Portable radio was detected