



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 ANNUAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 6492 Inspection Location: TOXL

A. Pre-Inspection

1. Items with Instrument:
  - Gas Cylinder Yes or No (If Yes, Lot # 14323080A1 Cyl. # 17)
  - Keys Yes or No
  - Power Cord Yes or No
2.  Download Data
3.  Upload Operator File
4.  Current Location Code: FYPD
5.  Battery Check
  - Was the external battery pack replaced? Yes or No
6.  O-Rings
  - Replaced Simulator O-Ring Yes or No
  - Replaced Breath Tube O-Ring Yes or No

B. General Setup and Checks:

1.  Diagnostics passed and instrument in "Ready" mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
  - Time Zone: CST or MDT (Time on test records will be in time zone circled).
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G).
  - Display: 1000 psi Regulator: 1000 psi
  - Display and Regulator ± 50 psi of each other Yes or No
  - Gas tank tare necessary? Yes or No
  - If Yes, display readings after tare (Level 3,M,C,G):
  - Display: \_\_\_\_\_ psi Regulator: \_\_\_\_\_ psi

C. Tests (Sign and attach test records):

1.  Configure simulator for the following test (Level 1,S).
  - Wet Calibration Check - Low AC (Level 1,C)
  - Known Value ≤ 0.03 AC: 0.020 AC
  - Sim. Ser #: MP3002
  - Lot #: 202301E
  - Exp. Date: 1/1/25
  - Results ± 0.005 of known AC

2.  Configure simulator for the following test (Level 1,S).  
 Wet Calibration Check - High AC (Level 1,C)  
 Known Value  $\geq$  0.25 AC: 0.400 AC  
 Sim. Ser #: MP6036  
 Lot #: 202212C  
 Exp. Date: 08Dec24  
 Results  $\pm$  5% AC of known AC
3.  Configure dry gas standard for the remaining tests (Level 1,S).  
 Known Value : 0.080 AC  
 Gas Cylinder Lot #: 01923080A3  
 Cylinder #: 18  
 Exp. Date: 215125
4.  Interferent Check (Level 1,B)  
 Known Value: 0.10 AC + 0.05% Acetone  
 Sim. Ser #: DR7344  
 Lot #: ICS 8  
 Exp. Date: N/A  
 Display reads "Interferent Detect"
5.  RFI Check (CMS Mode)  
 Display reads "RFI Detect"
6.  Dry Calibration Check (Level 1,C)  

|        |              |         |              |              |              |
|--------|--------------|---------|--------------|--------------|--------------|
| Test 1 | <u>0.081</u> | Test 4  | <u>0.081</u> | Test 7       | <u>0.081</u> |
| Test 2 | <u>0.081</u> | Test 5  | <u>0.081</u> | Test 8       | <u>0.081</u> |
| Test 3 | <u>0.081</u> | Test 6  | <u>0.081</u> | Test 9       | <u>0.081</u> |
|        |              | Average |              | <u>0.081</u> |              |

 Results  $\pm$  0.005 AC of known AC

D. Remarks/Maintenance: N/A

---



---

Instrument is acceptable to be used in the field.  Yes or No  
 If No, state reason(s) why: \_\_\_\_\_

If Yes, change location code back to A.4.

[Signature]  
 Inspector Signature

10May2024  
 Date

[Signature]  
 Reviewer

10May2024  
 Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-006492  
Location = TOXL                8164.14.00 09/16  
05/10/2024                      08:36

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$^&\*()\_+?

abcdefghijklmnopqrstuvwxy1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?         No  
Display Volume?             No  
Disable On Memfull?         Yes  
# of Print Copies?          1  
Select Std (D/W/I)?         Dry  
Standard Value?             0.080  
Standard Lot #?             14323080A1  
Standard Cyl #?             17  
Standard Expiration?        06/05/2025  
Oper No?                    133237

Flow Cal. Date:             08/18/2015  
Slope                        669  
Intercept                   -631515

IR Calibration Date:        08/18/2015  
                              3um                    9um

|                    |        |        |
|--------------------|--------|--------|
| 0th Coef(*100):    | -41685 | -23145 |
| 1st Coef(*100):    | 269226 | 142194 |
| 2nd Coef(*100):    | 2671   | 1239   |
| H2O adj(mg/l*10k): | 872    | 437    |

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ANNA NAREHOOD

Remarks: *print Test*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006492  
Location = TOXL      8164.14.00 09/16  
05/10/2024      09:13

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 09:14 |
| 02 Std. Sol. | 0.022 | 09:15 |
| 03 Room Air  | 0.000 | 09:16 |
| 04 Std. Sol. | 0.021 | 09:16 |
| 05 Room Air  | 0.000 | 09:17 |
| 06 Std. Sol. | 0.021 | 09:18 |
| 07 Room Air  | 0.000 | 09:18 |

08 Sim Temp = 34.0°C

Simul Ser No = MP3002

Std Sol No = 202301E

County = 08

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *LOW AC Check- 0.020AE*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006492  
Location = TOXL      8164.14.00 09/16  
05/10/2024      09:19

WET CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 09:20 |
| 02 Std. Sol. | 0.397 | 09:20 |
| 03 Room Air  | 0.000 | 09:21 |
| 04 Std. Sol. | 0.400 | 09:22 |
| 05 Room Air  | 0.000 | 09:22 |
| 06 Std. Sol. | 0.400 | 09:23 |
| 07 Room Air  | 0.000 | 09:23 |


08 Sim Temp = 34.0°C

Simul Ser No = MP6036

Std Sol No = 202212C

County = 08

Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: ~~Inte~~ High Ac Check - 0.400AE

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-006492  
Location = TOXL                              8164.14.00 09/16  
05/10/2024    09:26

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Room Air      | 0.000 | 09:27 |
| 02 *Subject Test | INT*  | 09:27 |
| 03 Room Air      | 0.000 | 09:28 |

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                                      Weight = NA  
Test = OTH                                      Cit = INTERFERENT CK  
Dr. Lic. = ND/TES989643  
Lot No = 01923080A3  
Cyl No = 18  
Expiration Date = 02/05/2025  
County = 08                                      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Interferent check*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-006492  
Location = TOXL                                      8164.14.00 09/16  
05/10/2024    09:29

| Test             | AC    | Time  |
|------------------|-------|-------|
| 01 Diagnostic    | OK    | 09:29 |
| 02 Room Air      | 0.000 | 09:30 |
| 03 *Subject Test | RFI*  | 09:30 |
| 04 Room Air      | 0.000 | 09:30 |

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                                      Weight = NA  
Test = OTH                                      Cit = RFI CHECK  
Dr. Lic. = ND/TES989643  
Lot No = 01923080A3  
Cyl No = 18  
Expiration Date = 02/05/2025  
County = 08                                      Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA NAREHOOD

Remarks: RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006492  
Location = TOXL      8164.14.00 09/16  
05/10/2024      09:31

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:31 |
| 02 Std. Gas | 0.081 | 09:32 |
| 03 Room Air | 0.000 | 09:32 |
| 04 Std. Gas | 0.081 | 09:33 |
| 05 Room Air | 0.000 | 09:33 |
| 06 Std. Gas | 0.081 | 09:33 |
| 07 Room Air | 0.000 | 09:34 |

Lot No = 01923080A3  
Cyl No = 18  
Exp Date = 02/05/2025  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Dry cal check #1-3*

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006492  
Location = TOXL      8164.14.00 09/16  
05/10/2024      09:34

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:35 |
| 02 Std. Gas | 0.081 | 09:35 |
| 03 Room Air | 0.000 | 09:36 |
| 04 Std. Gas | 0.081 | 09:36 |
| 05 Room Air | 0.000 | 09:37 |
| 06 Std. Gas | 0.081 | 09:37 |
| 07 Room Air | 0.000 | 09:38 |

Lot No = 01923080A3  
Cyl No = 18  
Exp Date = 02/05/2025  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: Dry cal check  
#4-6

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006492  
Location = TOXL      8164.14.00 09/16  
05/10/2024      09:38

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 09:39 |
| 02 Std. Gas | 0.081 | 09:39 |
| 03 Room Air | 0.000 | 09:40 |
| 04 Std. Gas | 0.081 | 09:40 |
| 05 Room Air | 0.000 | 09:41 |
| 06 Std. Gas | 0.081 | 09:41 |
| 07 Room Air | 0.000 | 09:42 |

Lot No = 01923080A3  
Cyl No = 18  
Exp Date = 02/05/2025  
County = 08      Oper No. = 133237



Operator Signature  
ANNA NAREHOOD

Remarks: *Dry cal check #7-9*

Form 106-I8000