

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-00 66 79	Instrument Location HILL				
Reason for Install/Repair					
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change					
(Specify) annual Inspection at Crime Laboratory					
THE PROPERTY OF THE PROPERTY O	<u></u>				
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
1 6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).					
7. Repair and/or Maintenance Performed (if any):					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form				
Field Inspector Signature	Date 5 / 10 / 2024				
Crime Laboratory Use Only					
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct county.					
Reviewed/Certified By	Certified Date 13May2024				

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006679
Location = HILL 8164.14.00 09/16
05/10/2024 20:10

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes Yes Display Third Digit? Inhib Printer(Y/N)? No Display Volume? NoDisable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 14323080A1 Standard Expiration? 06/05/2025 131630 Oper No?

Flow Cal. Date: 09/06/2016 Slope 701 Intercept -682926

IR Calibration Date: 09/06/2016
3um 9um

 0th Coef(\*100):
 -18983
 -23165

 1st Coef(\*100):
 280986
 134481

 2nd Coef(\*100):
 2443
 1290

 H2O adj(mg/l\*10k):
 569
 553

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature
AUSTIN BARIL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HILL
05/10/2024

Alcohol Analyzer SN 80-006679 8164.14.00 09/16 20:10

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:11
02 Std. Gas	0.081	20:12
03 Room Air	0.000	20:12
04 Std. Gas	0.081	20:12
05 Room Air	0.000	20:13
06 Std. Gas	0.081	20:13
07 Room Air	0.000	20:14

Lot No = 14323080A1

Cyl No = 67

Exp Date = 06/05/2025

County = 49

Oper No. = 131630

Operator Signature AUSTIN BARIL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006679 Location = HILL 8164.14.00 09/16 05/10/2024

20:16

Test			AC	Time
01	Room	Air	RFI*	20:20
02	Room	Air	0.000	20:20

\*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI TEST

Sub DOB = 01/01/2001

Sub Sex = Unknown

Weight = NA

Test = OTH

Cit = NA

Dr. Lic. = ND/NA

Lot No = 14323080A1

Cyl No = 67

Expiration Date = 06/05/2025

County = 49

Oper No. = 131630

Operator Signature AUSTIN BARIL

Remarks:

Form 106-I8000