

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59281 (06/2018)

Serial Number	80-004/89	Instrument Location NTPD
Reason for Ins	tall/Repair	
☐ Install	After Receiving From Crime Laboratory Install	I After Location Change
	(Specify) annual Inspection at Crim	ac Laboratory
Check When D	Done:	
	rge Protector Installed/Property Grounded.	
	lephone Line Connected to Intoxilyzer® 8000.	
	eath Tube Heated.	
	ter Preliminary Data (i.e. Date, Time, DST (Y), and Lo	ocation; Level 2, Function E).
∑ 5. Sc	an/Enter Gas Cylinder Information (Level 1, Function	n S).
	n Tests:	
	Print Test (Level1, Function P).	
	ACA Test (Level 1, Function C).	
	Radio Frequency Interference (RFI) Test (CMS Moo	de or Level 1, Function B or C; Key Radio During Test).
	epair and/or Maintenance Performed (if any):	me_
	omplete the Top Portion of the Intoxilyzer® Record (S	SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
	e Previous Intoxilyzer® Record (SFN504096, Form 1	20-G) at the Intoxilyzer® Location at the Agency.
10. Se	end the Following to the Crime Laboratory: Complete	ed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
10	04-G), Print Test, ACA Test, and RFI Test.	
Field Inspecto	r Signature	Date
Ch	- (3/	5/17/24
Crime Laborat	tony Use Only	
	the instrument is approved	d to be used for the analysis of breath to determine alcohol concentration fro
the date the F	ield Inspector performed the installation. This record	of the at the office of Attorney deficient, crime
County of Bur	leigh, North Dakota, is certified to be a true and corre	ect copy of the documents received.
Reviewed/Cer	rtified By	Certified Date 20May2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004189
Location = NTPD	8164.14.00 09/16
05/03/2024	16:29

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer (Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	14323080A1
Standard Cyl #?	72
Standard Expiration?	06/05/2025
Oper No?	133123
Flow Cal. Date: Slope	02/02/2011 628
<u>F</u>	020

Flow Cal. Date:	02/02/2011
Slope	628
Intercept	-452608

IR Ca	libration	Date	:	05/19/20	16
			3um		9um
	oef(*100):		-17767		-18761
	oef(*100):		263439		133606
2nd C	oef(*100):]	L960		1140
H20 a	di (ma/1*10	k): 3	341		362

****** Printer Test End *******

Operator Signature ORIE GOUDGE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = NTPD 05/03/2024

SN 80-004189 8164.14.00 09/16 16:46

DRY CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	16:46
02	Std.	Gas	0.079	16:47
03	Room	Air	0.000	16:47
04	Std.	Gas	0.079	16:48
05	Room	Air	0.000	16:48
06	Std.	Gas	0.080	16:48
07	Room	Air	0.000	16:49

Lot No = 14323080A1

Cyl No = 72

Exp Date = 06/05/2025

County = 31

Oper No. = 133123

Operator Signature ORIE GOUDGE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = NTPD 05/17/2024

SN 80-004189 8164.14.00 09/16 16:16

Test	AC	Time
01 Room Air	RFI*	
02 Room Air	0.000	16:18

*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI TEST

Sub DOB = 01/21/1975

Weight = NA Sub Sex = Unknown

Test = OTH

Cit = NA

Dr. Lic. = ND/NA Lot No = 14323080A1

Cyl No = 72

Expiration Date = 06/05/2025

County = 60

Oper No. = 133123

Operator Signature ORIE GOUDGE

Remarks:

Form 106-I8000