



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <b>80-006689</b>	Instrument Location <b>STARKSD/SINMCC DKPD</b>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u>Annual Inspection at Crime Laboratory</u>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <b>[Signature] #288</b>	Date <b>5/20/24</b>
--	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <b>[Signature]</b>	Certified Date <b>21May2024</b>
---	------------------------------------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006689  
Location = DKPD                    8164.14.00 09/16  
05/20/2024                            15:55

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer(Y/N)?                No  
Display Volume?                    No  
Disable On Memfull?                Yes  
# of Print Copies?                  1  
Select Std (D/W/I)?                Dry  
Standard Value?                    0.080  
Standard Lot #?                    14323080A4  
Standard Cyl #?                    40  
Standard Expiration?                06/05/2025  
Oper No?                            130396

Flow Cal. Date:                    05/29/2020  
Slope                                643  
Intercept                            -576193

IR Calibration Date:                05/29/2020  
    3um                                9um

	3um	9um
0th Coef(*100):	-25870	-28277
1st Coef(*100):	287704	139899
2nd Coef(*100):	26	744
H2O adj(mg/l*10k):	293	328

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ERIN QUINN

Remarks:

Form 106-I8000

The Undersigned, having custody of the original public record, certifies the information hereon is a true and correct copy of the original document maintained as a part of the activities of this office.

Signature/Title Erin Quinn - SOT  
Office/Agency NDATP


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-006689  
Location = DKPD                              8164.14.00 09/16  
05/20/2024    15:55

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:56
02 Std. Gas	0.079	15:56
03 Room Air	0.000	15:57
04 Std. Gas	0.079	15:57
05 Room Air	0.000	15:58
06 Std. Gas	0.078	15:58
07 Room Air	0.000	15:59

Lot No = 14323080A4  
Cyl No = 40  
Exp Date = 06/05/2025  
County = 45    Oper No. = 130396



Operator Signature  
ERIN QUINN

Remarks:

Form 106-I8000

The Undersigned, having custody of the original public record, certifies the information hereon is a true and correct copy of the original document maintained as a part of the activities of this office.

Signature/Title Erin Quinn SGT  
Office/Agency NDHP

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006689  
Location = DKPD      8164.14.00 09/16  
05/20/2024      16:00

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	16:00
02 Room Air	0.000	16:00

\*Invalid Test  
Inhibited - RFI

Lot No = 14323080A4  
Cyl No = 40  
Exp Date = 06/05/2025  
County = 45      Oper No. = 130396



Operator Signature  
ERIN QUINN

Remarks:

Form 106-I8000

The Undersigned, having custody of the original public record, certifies the information hereon is a true and correct copy of the original document maintained as a part of the activities of this office.

Signature/Title Erin Quinn SGT  
Office/Agency NDHP