

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 70 95 Instrument Location Calves PC)				
Reason for Install/Repair					
☐ Install After Receiving From Crime Laboratory ☐ Install After Location Change					
* Other (Specify) annual Inspection at Crime Laboratory					
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio D	uring Test).				
7. Repair and/or Maintenance Performed (if any):					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the	Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Ag	ency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.					
Field Inspector Signature Makhau C. Rama	Date 5/22/24				
Crime Laboratory Use Only	/ /				
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.					
Reviewed/Certified By	Certified Date				

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = OKPD 8164.16.00 09/18
05/22/2024 15:31

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 14323080A1 63 Standard Expiration? 06/05/2025 130690 Oper No?

Flow Cal. Date: 09/27/2018 Slope 683 Intercept -637452

IR Calibration Date: 06/21/2023

3um 9um

Oth Coef(*100): -38728 -15234

1st Coef(*100): 250485 139595

2nd Coef(*100): 2553 1269

H2O adj(mg/l*10k): 921 383

****** Printer Test End *******

Operator Signature MATTHEW O'BRIEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-007095 Location = OKPD 8164.16.00 09/18 05/22/2024

15:32

DRY CAL CHECK

Test			AC	Time	
	01	Room	Air	0.000	15:33
	02	Std.	Gas	0.080	15:34
	03	Room	Air	0.000	15:34
	04	Std.	Gas	0.080	15:35
	05	Room	Air	0.000	15:35
	06	Std.	Gas	0.080	15:36
	07	Room	Air	0.000	15:36

Lot No = 14323080A1

Cyl No = 63

Exp Date = 06/05/2025

County = 11

Oper No. = 130690

Operator Signature MATTHEW O'BRIEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = OKPD 05/22/2024

SN 80-007095 8164.16.00 09/18 15:37

Test			AC	Time
01	Room	Air	RFI*	15:39
02	Room	Air	0.000	15:40

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 01/27/1991

Sub Sex = Male

Weight = NA

Test = OTH

Cit = NA

Dr. Lic. = ND/NA

Lot No = 14323080A1

Cyl No = 63

Expiration Date = 06/05/2025

County = 11

Oper No. = 130690

Operator Signature MATTHEW O'BRIEN

Remarks:

Form 106-I8000