

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	1 /
80-005943	NCSO- NE Ville Hosp	ia l
Reason for Install/Repair		
Install / titol / teoring	er Location Change	
Other (Specify) Annual Inspection at	Onme Laboratory	
	U	
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locat	tion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
€ 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode of	r Level 1, Function B or C; Key Radio Dur	ring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFNs	50496, Form 120-G) and Place it by the In	toxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-	G) at the Intoxilyzer® Location at the Age	ncy.
10. Send the Following to the Crime Laboratory: Completed In	toxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
104-G), Print Test, ACA Test, and RFI Test.		
Field Inspector Signature		Date
Madela		5/24/2024
Crime Laboratory Use Only This installation has been reviewed and the instrument is approved to	he used for the analysis of breath to deter	rmine alcohol concentration from
the date the Field Inspector performed the installation. This record on	file at the Office of Attorney Ocherat, Other	ne Laboratory Division, in the
County of Burleigh, North Dakota, is certified to be a true and correct of	copy of the documents received.	
Reviewed/Certified By		Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005943
Location = NCSO 8164.14.00 09/16
05/24/2024 11:05

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 28423080A3 Standard Cyl #? 42 Standard Expiration? 11/05/2025 Oper No? 130917

Flow Cal. Date: 06/20/2013 Slope 686 Intercept -735246

IR Calibration Date: 08/26/2016 3um 9um

Oth Coef(*100): -13589 -19547 1st Coef(*100): 284566 139614 2nd Coef(*100): 2717 1336 H2O adj(mg/l*10k): 381 404

****** Printer Test End *******

Operator Signature
JESSE MADCHE

Remarks:

Print test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NCSO
05/24/2024

Alcohol Analyzer SN 80-005943 8164.14.00 09/16 11:05

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:06
02 Std. Gas	0.080	11:06
03 Room Air	0.000	11:07
04 Std. Gas	0.080	11:07
05 Room Air	0.000	11:08
06 Std. Gas	0.080	11:08
07 Room Air	0.000	11:09

Lot No = 28423080A3

Cyl No = 42

Exp Date = 11/05/2025

County = 32

Oper No. = 130917

Operator Signature JESSE MADCHE

Remarks:

Cal check

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NCSO
05/24/2024

Alcohol Analyzer SN 80-005943 8164.14.00 09/16 11:11

DRY CAL CHECK

Te	est ———	AC	Time
	Room Room	 RFI* 0.000	11:12

*Invalid Test Inhibited - RFI

Lot No = 28423080A3

Cyl No = 42

Exp Date = 11/05/2025

County = 32

Oper No. = 130917

Operator Signature JESSE MADCHE

Remarks:

Form 106-I8000

RFI test