

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	2. /		
80-004953	HIPTC/C	Lugby		
Reason for Install/Repair		0 /		
	er Location Change			
& Other (Specify) Annual Inspection at Crime Laboratory				
Check When Done:	0			
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
2 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.				
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.				
Field Inspector Signature		Date		
- Charles		6-3-24		

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By Certified	
an the	tified Date 04June2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004953 Location = RUGB 8164.14.00 09/16 06/03/2024 17:01 abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 14323080A1 Standard Cyl #? 69 Standard Expiration? 06/05/2025 Oper No? 130683 Flow Cal. Date: 04/20/2018 Slope 609 Intercept -462839 IR Calibration Date: 04/20/2018 3um 9um Oth Coef(*100): -22916 -27195 1st Coef(*100): 261768 2nd Coef(*100): 1749 134580 1216 540 H2O adj(mg/l*10k): 496

********** Printer Test End **********

Tast

Operator Signature SHAINE TOWE

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004953
Location = RUGB	8164.14.00 09/16
06/03/2024	17:01

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.080 0.000 0.081 0.000	17:02 17:02 17:03 17:03 17:04 17:04 17:04

The

Operator Signature SHAINE TOWE

Remarks:

Form 106-18000

ę.

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004953Location = RUGB8164.14.00 09/1606/03/202417:05

DRY CAL	CHECK	
Test	AC	Time
01 Room Air 02 Room Air	RFI* RFI*	17:06 17:06
*Invalid Test Inhibited - RFI		
Lot No = 14323080A1 Cyl No = 69 Exp Date = 06/05/2025 County = 35	Oper No.	= 130683

Operator Signature SHAINE TOWE

Remarks:

Form 106-18000