

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 20-006681	Instrument Location	25
Reason for Install/Repair		
☐ Install After Receiving From Crime Laboratory ☐ Install After	r Location Change	
Other (Specify)		
· .		
Check When Done:		
Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):	he	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G)) at the Intoxilyzer® Location at the Age	ncy.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	xilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
	Ψ ,	
Field Inspector Signature		Date 6-4-24
	M.	
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on file County of Burleigh, North Dakota, is certified to be a true and correct coperation.	e at the Office of Attorney General, Crin	
Reviewed/Certified By		Certified Date
/// / / / / / / / / / / / / / / / / / /		05 June 2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006681
Location = CARS	8164.14.00 09/16
06/04/2024	18:43

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	14323080A1
Standard Cyl #?	65
Standard Expiration?	06/05/2025
Oper No?	132318

Flow Cal.	Date:	05/27/2020
Slope		651
Intercept		-622365

IR Calibration Date	e: 05/27/2	020
	3um	9um
0th Coef(*100):	-10420	-19573
1st Coef(*100):	269401	137403
2nd Coef(*100):	1095	923
H2O adj(mg/l*10k):	127	287

****** Printer Test End *******

Operator Signature TERRY DILLON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = CARS
06/04/2024

Alcohol Analyzer SN 80-006681 8164.14.00 09/16 18:44

DRY CAL CHECK

Test		AC	Time
01 Room	Air	0.000	18:45
02 Std.	Gas	0.078	18:45
03 Room	Air	0.000	18:45
04 Std.	Gas	0.079	18:46
05 Room	Air	0.000	18:46
06 Std.	Gas	0.078	18:47
07 Room	Air	0.000	18:47

Lot No = 14323080A1

Cyl No = 65

Exp Date = 06/05/2025

County = 19

Oper No. = 132318

Operator Signature TERRY DILLON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006681
Location = CARS	8164.14.00 09/16
06/04/2024	18:51

Test	AC	Time
01 Room Air	0.000	18:52
02 *Subject Test	RFI*	18:53
03 Room Air	RFI*	18:53

*Invalid Test Inhibited - RFI

Sub Name = DILLON, TERRY LEE

Sub DOB = 05/14/1959

Sub Sex = Male Weight = NA Test = DUI Cit = TEST

Dr. Lic. = ND/DIL598214

Lot No = 14323080A1

Cyl No = 65

Expiration Date = 06/05/2025

County = 19 Oper No. = 132318

Operator Signature
TERRY DILLON

Remarks:

Form 106-I8000