

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 003067	Instrument Location Grafton	Pd
Reason for Install/Repair		
☐ Install After Receiving From Crime Laboratory ☐ Install After	er Location Change	
M Other (Specify) annual Inspection at	Crime Laboratory	
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
e 6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50	0496, Form 120-G) and Place it by the li	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G	s) at the Intoxilyzer® Location at the Age	ncy.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature		Date (0.6.24
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be	a used for the analysis of broath to dotor	mine alcohol concentration from
the date the Field Inspector performed the installation. This record on fill County of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorney General, Crin	
Reviewed/Certified By	0	Certified Date
(/Me / Inle	X	06June2024

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = GFTN 8164.14.00 09/16
06/06/2024 13:27

********* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01923080A3 Standard Cyl #? 16 02/05/2025 Standard Expiration? Oper No? 020305

Flow Cal. Date: 04/17/2020 Slope 609 Intercept -236750

IR Calibration Date: 04/17/2020
3um 9um

Oth Coef(*100): -19307 -18350
1st Coef(*100): 282277 142540
2nd Coef(*100): 1257 744
H2O adj(mg/l*10k): 302 214

******* Printer Test End *******

Operator Signature
ANTHONY DUMAS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFTN
06/06/2024

Alcohol Analyzer SN 80-003067 8164.14.00 09/16 13:31

DRY CAL CHECK

Test		AC	Time
01 Room F	Air 0	.000	13:32
02 Std. 0	Gas 0	.080	13:32
03 Room A	Air 0	.000	13:33
04 Std. 0	Gas 0	.079	13:33
05 Room A	Air 0	.000	13:34
06 Std. 0	as 0	.079	13:34
07 Room A	Air 0	.000	13:34

Lot No = 01923080A3

Cyl No = 16

Exp Date = 02/05/2025

County = 50

Oper No. = 020305

Operator Signature ANTHONY DUMAS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFTN
06/06/2024

Alcohol Analyzer SN 80-003067 8164.14.00 09/16 13:35

$T\epsilon$	est		AC Tim	ne —
01	Room	Air	RFI* 13:3	6
02	Room	Air	RFI* 13:3	6

*Invalid Test Inhibited - RFI

Sub Name = DUMAS, ANTHONY GERALD

Sub DOB = 01/07/1983

Sub Sex = Male

Weight = 230

Cit = 1435

Test = OTH

Dr. Lic. = ND/DUM830566

Lot No = 01923080A3

Cyl No = 16

Expiration Date = 02/05/2025

County = 50

Oper No. = 020305

Operator Signature ANTHONY DUMAS

Remarks:

Form 106-I8000