

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 5 3 5 7	Instrument Location	
Reason for Install/Repair		
☐ Install After Receiving From Crime Laboratory ☐ Ins	stall After Location Change	
	at crime Laboratory	
Check When Done:	_	
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and	d Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Fund	tion S).	
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
EN C. Radio Frequency Interference (RFI) Test (CMS	Mode or Level 1, Function B or C; Key Ra	dio During Test).
7. Repair and/or Maintenance Performed (if any):		
☑ 8. Complete the Top Portion of the Intoxilyzer® Record		
2 9. File Previous Intoxilyzer® Record (SFN504096, For	rm 120-G) at the Intoxilyzer® Location at t	he Agency.
10. Send the Following to the Crime Laboratory: Comp 104-G), Print Test, ACA Test, and RFI Test.	oleted Intoxilyzer® 8000 Installation and R	epair Checkout (SFN59281, Form
Field Inspector Signature		Date 6/21/2024
Control of the Control		
Crime Laboratory Use Only This installation has been reviewed and the instrument is appropriate to the control of the control	oved to be used for the analysis of breath	to determine alcohol concentration from
the date the Field Inspector performed the installation. This re County of Burleigh, North Dakota, is certified to be a true and county of Burleigh.	cord on the strife office of Adolle's ocue	ral, Crime Laboratory Division, in the
Reviewed/Certified By	nl	Certified Date 24June2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005357
Location = GRAF	8164.14.00 09/16
06/20/2024	17:14

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()__+?

Current Instru	ment Setup
----------------	------------

Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	01923080A3
Standard Cyl #?	21
Standard Expiration?	02/05/2025
Oper No?	132597
_	

Flow Cal.	Date:	05/16/2018
Slope		573
Intercept		140460

IR Calibration Date	э: (06/21/2023
	3um	9um
		
0th Coef(*100):	-6039	-13097
1st Coef(*100):	257461	136071
2nd Coef(*100):	3187	1488
H20 adj(mg/l*10k):	315	418

****** Printer Test End *******

Operator Signature JACOB NIKKO

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = GRAF 8164.14.00 09/16 06/20/2024

SN 80-005357 17:16

DRY CAL CHECK

Test	AÇ	Time
01 Room Air	0.000	17:16
02 Std. Gas	0.080	17:17
03 Room Air	0.000	17:17
04 Std. Ģas	0.080	17:18
05 Room Air	0.000	17:18
06 Std. Gas	0.080	17:19
07 Room Air	0.000	17:19

Lot No = 01923080A3

Cyl No = 21

Exp Date = 02/05/2025

County = 50

Oper No. = 132597

Operator Signature JACOB NIKKO

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GRAF
06/20/2024

Alcohol Analyzer SN 80-005357 8164.14.00 09/16 17:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:21
02 Reference	RFI*	17:21
03 Room Air	0.000	17:21

*Invalid Test Inhibited - RFI

Lot No = 01923080A3

Cyl No = 21

Exp Date = 02/05/2025

County = 50

Oper No. = 132597

Operator Signature JACOB NIKKO

Remarks:

Form 106-I8000