

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

	INTOXILYZER® 8000 CALIBRATION	N
Intoxilyzer [®] 8	000 Serial Number: 80-00 5362 Calibratic	on Location: TDXL
A. Pre-Ca 1. 2. 3. 4.	alibration Download Data Upload Operator File Battery Check Was the external battery pack replaced? O-Rings Replaced Simulator O-Ring Replaced Breath Tube O-Ring	Yes or No Yes or No Yes or No
 B. Gener 1. 2. 3. 4. 5. 	 Setup and Checks: Diagnostics passed and instrument in "Ready" model Breath tube heated Date, time and location code (Level 2,E). Re-set is Time Zone: CST or MDT (Time on test record) Print test (Level 1,P). Sign and attach test record. Print test (Level 1,P). Sign and attach test record. Tank monitor (Level 3,D,G). Display: <u>100</u> psi Regulator: <u>150</u> Display and Regulator ± 50 psi of each other Gas tank tare necessary? If Yes, display readings after tare (Level 3,M,G Display: psi Regulator: 	if necessary. ds will be in time zone circled). yes or No Yes or No C,G):
C. Tests 1.	(Sign and attach test records): Configure simulator for the following test (Level 1 Wet Calibration Check - Low AC (Level 1,C) Known Value ≤ 0.03 AC: <u>0.020</u> AC Sim. Ser #: <u>MP300A</u> Lot #: <u>203301E</u> Exp. Date: <u>Dunas</u>	,S).

Laboratory Unit: Toxicology Unit - Breath Alcohol Section Approved By: Laboratory Director UNCONTROLLED WHEN PRINTED

Intoxilyzer 8000 Calibration

Results ± 0.005 of known AC

2.	Configure simulator for the following test (Level 1,S). Wet Calibration Check - High AC (Level 1,C) Known Value ≥ 0.25 AC: 0.400 AC Sim. Ser #: <u>NP6036</u> Lot #: <u>202212C</u> Exp. Date: <u>08 Dec24</u> M Results ± 5% AC of known AC
3.	Configure dry gas standard for the remaining tests (Level 1,S).
	Known Value : 0.080 AC
	Gas Cylinder Lot #: <u> </u>
	Cylinder #: 41 Exp. Date: 11 5 25
4.	Interferent Check (Level 1,B)
	Known Value: 0.10 AC + 0.05% Acetone Sim. Ser #:
	Lot #: ICS 8
	Exp. Date:NA
	Display reads "Interferent Detect"
5	RFI Check (CMS Mode)
	Display reads "RFI Detect"
6	Dry Calibration Check (Level 1,C)
0.	Test 1 0,081 Test 4 0.081 Test 7 0.082
	Test 1 0.081 Test 4 0.081 Test 7 0.082 Test 2 0.082 Test 5 0.081 Test 8 0.081
	Test 3 0.082 Test 6 0.081 Test 9 0.081
	Average 0.08
	Results ± 0.005 AC of known AC
-	NIO
Rema	arks/Maintenance: <u>NH</u>
	4. 5. 6.

Instrument is acceptable to be used in the field. (Yes)or No

If No, state reason(s) why:_

Calibration Analyst Signature

andle Postschiller Reviewer

Intoxilyzer 8000 Calibration Laboratory Unit: Toxicology Unit - Breath Alcohol Section Approved By: Laboratory Director UNCONTROLLED WHEN PRINTED

27 June 2024 Date

27 Sune rory Date

Document ID: 11871 Revision: 3 Status: Published Date Approved: 02/29/2024 Page 2 of 2

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005362
Location = TOXL	8164.14.00 09/16
06/27/2024	10:56

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

<pre>Inhib Printer(Y/N)? Display Volume? Disable On Memfull? # of Print Copies? Select Std (D/W/I)? Standard Value? Standard Lot #? Standard Cyl #? Standard Expiration? Oper No?</pre>	Enabled DABACABA Yes Yes No No Yes 1 Dry 0.080 14323080A4 15 06/05/2025 133237
	05/22/2020 659
Slope Intercept	-620990
IR Calibration Date:	05/26/2020
3um	9um
0th Coef(*100):-208831st Coef(*100):2655032nd Coef(*100):936	3 -25370

379

Operator Signature ANNA NAREHOOD

Remarks: PNAt Test

H2O adj(mg/l*10k): 391

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005362Location = TOXL8164.14.00 09/1606/27/202410:56

	T CAL CHECK		
Time	AC	est	Τe
10:57	0.000	Room Air	01
10:59	0.019	Std. Sol.	02
11:00	0.000	Room Air	03
11:00	0.019	Std. Sol.	04
11:01	0.000	Room Air	05
11:02	0.019	Std. Sol.	06
11:02	0.000	Room Air	07

08 Sim Temp = 34.0°C

Operator Signature ANNA NAREHOOD

Remarks: LOW AC Check - 0.020Ac

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005362Location = TOXL8164.14.00 09/1606/27/202411:03

		WET CAL CHECK	
Te	est	AC	Time
01	Room Air	0.000	11:03
02	Std. Sol.	0.379	11:04
03	Room Air	0.000	11:05
04	Std. Sol.	0.380	11:06
05	Room Air	0.000	11:06
06	Std. Sol.	0.381	11:07
07	Room Air	0.000	11:07
08	Sim Temp =	= 34.0°C	

Operator Signature ANNA NAREHOOD

Remarks:

High Accheek - 0.400Ac

CMI, Inc. Intoxily North Dakota Model Location = TOXL 06/27/2024	8000 SN	Analyzer 80-005362 .00 09/16 11:23
Test	AC	Time
01 Room Air	0.000	11:24
02 *Subject Test	INT*	11:25
03 Room Air	0.000	11:25

*Invalid Test Interferent Detected

Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998Sub Sex = Male Weight = NA Cit = INTERFERENT CK Test = OTHDr. Lic. = ND/TES989643 Lot No = 28423080A3Cyl No = 41Expiration Date = 11/05/2025County = 08Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: Interferent check

CMI, Inc. Intoxily North Dakota Model Location = TOXL 06/27/2024		Alcohol Analyzer SN 80-005362 8164.14.00 09/16 11:26
Test	AC	Time
01 Diagnostic 02 Room Air 03 Room Air	OK RFI 0.00	

*Invalid Test Inhibited - RFI

Sub Name = TEST, DONOR2 NONE Sub DOB = 07/25/1998 Sub Sex = Male Weight = NA Test = OTH Cit = RFI CHECK Dr. Lic. = ND/TES989643 Lot No = 28423080A3 Cyl No = 41 Expiration Date = 11/05/2025 County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature ANNA NAREHOOD

RFI Cherk

Remarks:

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005362
Location = TOXL	8164.14.00 09/16
06/27/2024	11:28

Г

	Ĺ	DRY CAL CHECK	
Test		AC	Time
Room	Air	0.000	11:28
Std.	Gas	0.081	11:28
Room	Air	0.000	11:29
Std.	Gas	0.082	11:29
Room	Air	0.000	11:30
Std.	Gas	0.082	11:30
Room	Air	0.000	11:31
	Room Std. Room Std. Room Std.	Room Air Std. Gas Room Air Std. Gas Room Air Std. Gas Room Air	ACRoom Air0.000Std. Gas0.081Room Air0.000Std. Gas0.082Room Air0.000Std. Gas0.082

Lot No = 28423080A3 Cyl No = 41 Exp Date = 11/05/2025 County = 08 Oper No. = 133237

Cal

Operator Signature ANNA NAREHOOD

Check #1-3

Remarks: DM

Form 106-18000

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005362Location = TOXL8164.14.00 09/1606/27/202411:31

	DRY CAL CHECK			
Time	AC		est	Τe
11:32	0.000	Air	Room	01
11:32	0.081	Gas	Std.	02
11:32	0.000	Air	Room	03
11:33	0.081	Gas	Std.	04
11:33	0.000	Air	Room	05
11:34	0.081	Gas	Std.	06
11:34	0.000	Air	Room	07

Lot No = 28423080A3 Cyl No = 41 Exp Date = 11/05/2025 County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

check #4-6

Remarks:

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005362
Location = TOXL	8164.14.00 09/16
06/27/2024	11:35

			DRY CAL (CHECK	
Test				AC	Time
01	Room	Air	(0.000	11:35
02	Std.	Gas	(0.082	11:36
03	Room	Air	(0.000	11:36
04	Std.	Gas	(0.081	11:37
05	Room	Air	(0.000	11:37
06	Std.	Gas	(0.081	11:38
07	Room	Air		0.000	11:38

Lot No = 28423080A3 Cyl No = 41 Exp Date = 11/05/2025 County = 08

Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: Dry Calcher #7-9