

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Instru	ment Location	
80-005958	MIFL	
Reason for Install/Repair		
Install After Receiving From Crime Laboratory Install After Loca	tion Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
☑ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Le	vel 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level	1, Function B or C; Key Radio Du	uring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496,	Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the	e Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyze 104-G), Print Test, ACA Test, and RFI Test.	r® 8000 Installation and Repair C	Checkout (SFN59281, Form
Field Inspector Signature		Date ,
with latter		07/18/2024
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be used	for the analysis of breath to doto	mine alcohol concentration from
the date the Field Inspector performed the installation. This record on file at the County of Burleigh, North Dakota, is certified to be a true and correct copy of	e Office of Attorney General, Cri	
Reviewed/Certified By		Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005958
Location = MIFL	8164.14.00 09/16
07/18/2024	16:39

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
<pre>Inhib Printer(Y/N)?</pre>	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	01524080A4
Standard Cyl #?	37
Standard Expiration?	03/05/2026
Oper No?	100913

Flow Cal. Date:	06/21/2013
Slope	716
Intercept	-760436

IR Calibration Date	e: 04/07/2	015
	3um	9um
0th Coef(*100):	-17179	-17694
1st Coef(*100):	259328	130890
2nd Coef(*100):	3473	1654
H20 adj(mq/l*10k):	425	422

****** Printer Test End *******

Operator Signature
JUSTIN COTE-KANNING

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MIFL
07/18/2024

Alcohol Analyzer SN 80-005958 8164.14.00 09/16 16:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:36
02 Std. Gas	0.079	16:36
03 Room Air	0.000	16:37
04 Std. Gas	0.078	16:37
05 Room Air	0.000	16:38
06 Std. Gas	0.078	16:38
07 Room Air	0.000	16:39

Lot No = 01524080A4

Cyl No = 37

Exp Date = 03/05/2026

County = 51

Oper No. = 100913

Operator Signature
JUSTIN COTE-KANNING

Remarks:

Form 106-I8000

intoxilyzen Test Record and Checklist NOOPS Crime Lab. Siv., Bishardk, NO 58501

emi, inc. Intexilyzer	Alconol Analyzer
Month Dakota Model 8000	SN 83-885958
Location = MIFL	8164.14.88 89/16
07/18/2024	16:41

Test	36	Tira
01 Room Air	RF:×	15:42
92 Room Air	RF:×	16:42

*invalid Test Inhibitea - RFl

Sub Name = TEST, TEST TEST Sub DOB = 01/01/1901

Sub Sex = Male

weight = 200 Cit = 0000000

Test = DUI

Or. Lic. = ND/TEST

Lot No = 0152408004

Eyi No = 37

Expiration Date = 93/05/2026

County = 51

Open No. = 100913

Operator Signature JUSTIN COTE-KANNING

Remarks:

Form 106-18000