

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-004199	Instrument Location Grafton	Po			
Reason for Install/Repair	9,1401.01				
Install After Receiving From Crime Laboratory 🔲 Install Aft	ter Location Change				
Other (Specify)					
Check When Done:					
	1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.					
3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode o	r Level 1, Function B or C; Key Radio Du	rring Test).			
7. Repair and/or Maintenance Performed (if any):					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	50496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0)	G) at the Intoxilyzer® Location at the Age	ency.			
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form			
Field Inspector Signature		Date 7.29.2024			
Crime Laboratory Use Only					
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on f County of Burleigh, North Dakota, is certified to be a true and correct co	file at the Office of Attorney General, Crir				
Reviewed/Certified By)	Certified Date 01Aug2024			

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004199
Location = GFTN 8164.14.00 09/16
07/29/2024 08:02

********* Printer Test *********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: DABACABA Start Test Sequence: Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? 14323080A4 Standard Lot #? 43 Standard Cyl #? 06/05/2025 Standard Expiration? 020305 Oper No?

Flow Cal. Date: 06/17/2019 Slope 593 Intercept -481610

IR Calibration Date: 04/18/2024 3um 9um

 0th Coef(*100):
 -25088
 -28296

 1st Coef(*100):
 269537
 131183

 2nd Coef(*100):
 2308
 1373

 H2O adj(mg/l*10k):
 357
 563

******* Printer Test End ********

Operator Signature
ANTHONY DUMAS

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFTN
07/29/2024

Alcohol Analyzer SN 80-004199 8164.14.00 09/16 08:09

DRY CAL CHECK

Test		AC	Time
01 Room A 02 Std. 0 03 Room A 04 Std. 0 05 Room A	Gas Air Gas Air	0.000 0.079 0.000 0.079 0.000	08:10 08:10 08:11 08:11
06 Std. 0		0.079 0.000	08:12 08:12

Lot No = 01923080A3

Cyl No = 16

Exp Date = 02/05/2025

County = 50

Oper No. = 020305

ANTHONY DUMAS

Remarks:

Form 106-I8000

Operator Signature

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFTN
07/29/2024

Alcohol Analyzer SN 80-004199 8164.14.00 09/16 08:24

Test		AC	Time
	Room Room	 RFI* RFI*	08:25 08:25

*Invalid Test Inhibited - RFI

Sub Name = GROSSBAUER, JOHN EDWARD

Sub DOB = 01/18/1993

Sub Sex = Male

Weight = 220

Test = HUI Cit = 12345667

Dr. Lic. = ND/GRO936482 Lot No = 01923080A3

Cyl No = 16

Expiration Date = 02/05/2025

County = 50

Oper No. = 020305

Operator Signature
ANTHONY DUMAS

Remarks:

Form 106-I8000