

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005957 Instrument Location
Reason for Install/Repair
Install After Receiving From Crime Laboratory Install After Location Change
Check When Done:
1. Surge Protector Installed/Property Grounded.
2. Telephone Line Connected to Intoxilyzer® 8000.
3. Breath Tube Heated.
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
6. Run Tests:
A. Print Test (Level1, Function P).
B. ACA Test (Level 1, Function C).
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any): UPDATE
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.
Field Inspector Signature Miketton Date 25/2024

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

<u> </u>	A)		
Reviewed/Certified By	- Mall	Anna Narehood	Certified Date 26Aug2024	

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005957 Location = MAND 8164.14.00 09/16 08/25/2024 02:59 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 14323080A1 Standard Cyl #? 73 Standard Expiration? 06/05/2025 Oper No? 079706 Flow Cal. Date: 06/21/2013 Slope 647 Intercept -613317 IR Calibration Date: 05/11/2015 3um 9um ------Oth Coef(*100):-166661st Coef(*100):2784782nd Coef(*100):3303 -18045 132949 1397 H2O adj(mg/l*10k): 379351 ********** Printer Test End **********

Operator Signature DANA MIDDLETON

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005957Location = MAND8164.14.00 09/1608/25/202403:06

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.077 0.000 0.077 0.000 0.077 0.000	03:06 03:07 03:07 03:08 03:08 03:09 03:09

Lot No = 14323080A1 Cyl No = 73 Exp Date = 06/05/2025 County = 30

Oper No. = 079706

Operator Signature DANA MIDDLETON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-005957 Location = MAND 8164.14.00 09/16 08/25/2024 03:11 Time Test AC 01 Room Air 0.000 03:13 02 *Subject Test RFI* 03:13 03 Room Air 03:14 0.000 *Invalid Test Inhibited - RFI Sub Name = RFI, RFI RFI Sub DOB = 11/11/1988Weight = 200Sub Sex = Male Test = DUI Cit = NA Dr. Lic. = ND/NALot No = 14323080A1Cyl No = 73Expiration Date = 06/05/2025Oper No. = 079706County = 30

Operator Signature DANA MIDDLETON

Remarks:

Form 106-18000