

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005939	Instrument Location Williston			
Reason for Install/Repair				
	er Location Change			
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
☑ 3. Breath Tube Heated.				
√ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	☑ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
☑ 6. Run Tests:				
A. Print Test (Level1, Function P).				
☑ B. ACA Test (Level 1, Function C).				
☑ C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).				
√ 7. Repair and/or Maintenance Performed (if any):   √ //-  //-  //-  //-  //-  //-  //-	<i>Y</i>			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.		
	6) at the Intoxilyzer® Location at the Age	ency.		
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form		
Field Inspector Signature		Date		
776		9/27/24		
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on ficulty of Burleigh, North Dakota, is certified to be a true and correct co	ile at the Office of Attorney General, Crir			
Reviewed/Certified By Anna Narehood	10.10	Certified Date		
Allia Naicilou / // // A	1 Chl. X	01Oct2024		

01Oct2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005939
Location = WILL	8164.14.00 09/16
09/27/2024	13:34

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01524080A4 Standard Cyl #? 37 Standard Expiration? 03/05/2026 Oper No? 130996 Flow Cal. Date: 06/18/2020

Flow Cal. Date: 06/18/2020 Slope 611 Intercept -592495

H20 adj (mg/1\*10k): 220

IR Calibration Date: 06/18/2020

3um 9um

Oth Coef(\*100): -15726 -19510

1st Coef(\*100): 296627 143165

2nd Coef(\*100): -369 494

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

240

Operator Signature PRESTON LANGER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WILL
09/27/2024

Alcohol Analyzer SN 80-005939 8164.14.00 09/16 13:35

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:36
02 Std. Gas	0.080	13:36
03 Room Air	0.000	13:37
04 Std. Gas	0.078	13:37
05 Room Air	0.000	13:38
06 Std. Gas	0.078	13:38
07 Room Air	0.000	13:39

Lot No = 01524080A4

Cyl No = 37

Exp Date = 03/05/2026

County = 53

Oper No. = 130996

Operator Signature PRESTON LANGER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005939
Location = WILL	8164.14.00 09/16
09/27/2024	13:40

Test			AC	Time
01	Room	Air	RFI*	13:41
02	Room	Air	RFI*	13:42

\*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI TEST

Sub DOB = 01/01/2001

Sub Sex = Unknown Weight = NA

Test = OTH

Cit = NA

Dr. Lic. = ND/ND

Lot No = 01524080A4

Cyl No = 37

Expiration Date = 03/05/2026

County = 53

Oper No. = 130996

Operator Signature PRESTON LANGER

Remarks:

Form 106-I8000