

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00 4934	Instrument Location	morten	Detention Conter		
Reason for Install/Repair	0				
Install After Receiving From Crime Laboratory					
☐ Other (Specify)					
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
X 3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
₩ 6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).					
7. Repair and/or Maintenance Performed (if any):					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) a	nd Place it by t	the Intoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.					
 10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test. 	toxilyzer® 8000 Instal	llation and Rep	air Checkout (SFN59281, Form		
Field Inspector Signature			Date		
May1. Wm			01/15/2025		
Crime Laboratory Use Only			1		
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct c	file at the Office of Att	torney General,	determine alcohol concentration from , Crime Laboratory Division, in the		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004936 8164.14.00 09/16 Location = BMDC 01/15/2025 14:16abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABA Display Prelim Rslt? Yes DABACABA Display Third Digit? Yes Inhib Printer(Y/N)? No No Display Volume? Yes Disable On Memfull? # of Print Copies? 1 # Of Princ copies.Select Std (D/W/I)?DryStandard Value?0.080 Standard Lot #? Standard Cyl #? 01524080A4 27 Standard Expiration? 03/05/2026 Oper No? 131856 01/18/2010 Flow Cal. Date: Slope 732 Intercept -767769 IR Calibration Date: 04/29/2024 3um 9um 0th Coef(*100): -24811 -25942 lst Coef(*100): 260347 2nd Coef(*100): 2582 130533 1344 H2O adj(mg/l*10k): 289 383

erator Signature

TAYLOR WRAY

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004936Location = BMDC8164.14.00 09/1601/15/202514:17

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas	0.000 0.079 0.000 0.080 0.000 0.081	14:17 14:18 14:18 14:19 14:19 14:19

Impl-h Operator Signature TAYLOR WRAY

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004936Location = BMDC8164.14.00 09/1601/15/202514:25

DRY CAL	CHECK	
Test	AC	Time
01 Room Air 02 Room Air	RFI* RFI*	14:25 14:25
*Invalid Test Inhibited - RFI		
Lot No = 01524080A4 Cyl No = 27 Exp Date = 03/05/2026 County = 08	Oper No.	= 131856

Operator Signature

TAYLOR WRAY

Remarks:

Form 106-18000