



**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION**

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-00 4945 Calibration Adjustment Location: TOXL

A. Pre-Adjustment

Replaced Simulator Return O-Ring Yes or (No)

B. Calibration Adjustment (Level 3,M,C,O)

1. Autocalibration Printout Attached

Max Power Res Value ≥ 10

Auto Range Res Value ≥ 4

2. Simulator Solutions for Calibration Adjustment

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000	NA-Milli-Q H ₂ O	NA-Milli-Q H ₂ O	MP3003
2	0.040	202303H	28 Mar 25	MP3061
3	0.080	202302B	14 Feb 25	MP3070
4	0.100	202304A	04 Apr 25	MP3057
5	0.300	202402C	14 Feb 26	MP3048

3. 0.080 AC Calibration Gas for H₂O Adjustment

Lot No. 01923080A3 Cyl No. 22 Exp. Date: 2/5/2025

4. Atmospheric Pressure

Displayed by Intoxilyzer® 8000 946 mbar

Adjusted to using barometer 956 mbar

Auto Calibration Report printout 956 mbar

Barometer Model 10510-922

Barometer Serial Number 230307250

Barometer Calibration Expiration Date 02 May 2025

5. Screen displayed "Calibration Success"

6. Calibration Adjustment Printout Attached

Solution 1 Avg % Abs ≤ 0.2500

Solution 2-5 REL STD DEV ≤ 3.000

Residual (g/210 L) values for solutions 1 - 5 ≤ 0.0020 for 3 μ m and 9 μ m channels

Dry Gas H₂O adjustment sum for 3 μm and 9 μm channels within ± 10
 3 μm 3348.3 (Ave.) + 461 (H₂O Adj.) = 3809.3
 9 μm 3444.0 (Ave.) + 365 (H₂O Adj.) = 3809.0

C. Is an Annual Inspection due for this instrument? Yes or No
 If Yes, complete Intoxilyzer 8000 Annual Inspection (Document ID: 11698)
 If No, complete Intoxilyzer 8000 Calibration (Document ID: 11871).

Remarks/Notes: Forgot to update date/time before beginning tests.

Amy Miller
 Analyst Signature

01/22/2025
 Date

Janelle Putschler
 Reviewer Signature

23 Jan 2025
 Date

~~EMEC~~ *TOXL*
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-004945
~~01/21/2025~~ *aen* ~~11/22/25~~ ~~10:15:41~~
01/22/2025 *11/22/25* *10:15:41*
Auto Calibration
Max Power Res Value = 26
Auto Range Res Value = 14

BMDC ~~TOX~~

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-004945

01/21/2025

22:32:41

01/22/2025

aww 1/22/25

10:15:41

Auto Calibration

pg 1 of 2

<<<<<			3um	>>>>>			<<<<<			9um	>>>>>		

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)			% Abs	(% Abs Ref)				
Sample #1	0.1360	(0.0010)		0.1950	(0.0010)								
Sample #2	0.0730	(0.0590)		0.1740	(0.0300)								
Sample #3	0.0990	(0.0590)		0.1940	(0.0160)								
Sample #4	0.1220	(0.0650)		0.2130	(0.0140)								
Avg % Abs	0.0980	(0.0610)		0.1937	(0.0200)								
STD DEV	0.0245	(0.0035)		0.0195	(0.0087)								
REL STD DEV	25.016	(5.679)		10.070	(43.589)								

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)								
Sample #1	0.8060	(0.0070)		1.5020	(0.0100)								
Sample #2	0.8300	(0.0220)		1.4940	(0.0470)								
Sample #3	0.8050	(0.0470)		1.4970	(0.0440)								
Sample #4	0.8150	(0.0370)		1.4960	(0.0400)								
Avg % Abs	0.8167	(0.0353)		1.4957	(0.0437)								
STD DEV	0.0126	(0.0126)		0.0015	(0.0035)								
REL STD DEV	1.541	(35.612)		0.102	(8.042)								

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)								
Sample #1	1.5410	(-0.0080)		2.7910	(-0.0090)								
Sample #2	1.4900	(0.0350)		2.7710	(0.0250)								
Sample #3	1.5220	(0.0350)		2.7990	(0.0280)								
Sample #4	1.5080	(0.0490)		2.7750	(0.0360)								
Avg % Abs	1.5067	(0.0397)		2.7817	(0.0297)								
STD DEV	0.0160	(0.0081)		0.0151	(0.0057)								
REL STD DEV	1.065	(20.377)		0.544	(19.167)								

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)								
Sample #1	1.8780	(-0.0080)		3.4540	(-0.0030)								
Sample #2	1.8740	(0.0280)		3.4320	(0.0310)								
Sample #3	1.8680	(0.0180)		3.4090	(0.0440)								
Sample #4	1.8810	(0.0330)		3.4530	(0.0440)								
Avg % Abs	1.8743	(0.0263)		3.4313	(0.0397)								
STD DEV	0.0065	(0.0076)		0.0220	(0.0075)								
REL STD DEV	0.347	(29.004)		0.641	(18.922)								

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1													
Sample	% Abs	(% Abs Ref)		% Abs	(% Abs Ref)								
Sample #1	5.2510	(-0.0060)		9.4540	(-0.0050)								
Sample #2	5.2400	(0.0620)		9.4060	(0.1080)								
Sample #3	5.2430	(0.0670)		9.4340	(0.1090)								
Sample #4	5.2370	(0.0940)		9.4290	(0.1450)								
Avg % Abs	5.2400	(0.0743)		9.4230	(0.1207)								
STD DEV	0.0030	(0.0172)		0.0149	(0.0211)								
REL STD DEV	0.057	(23.158)		0.158	(17.469)								

BMDC ^{TDXL}
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-004945
 01/21/2025 22:32:41
 01/22/2025 AEM 01/22/25 10:15:41
 Auto Calibration

<<<<< 3um >>>>>			<<<<< 9um >>>>>		
-----			-----		
Zero Order Coef	-260.09			-270.15	
First Order Coef	2642.94			1428.61	
Second Order Coef	25.36			12.31	
-----			-----		
Act	Fit	Residual	Act	Fit	Residual
(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)
0.000	-0.000	0.0000	0.000	0.000	-0.0001
0.040	0.040	-0.0002	0.040	0.040	0.0002
0.080	0.079	0.0006	0.080	0.080	0.0002
0.100	0.100	-0.0004	0.100	0.100	-0.0003
0.300	0.300	0.0000	0.300	0.300	0.0000
-----			-----		

<<<<< 3um >>>>>		<<<<< 9um >>>>>	
-----		-----	
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1			
Sample			
Sample #1	3355.00		3429.00
Sample #2	3375.00		3454.00
Sample #3	3375.00		3447.00
Sample #4	3295.00		3431.00
Avg	3348.3333		3444.0000
STD DEV	46.1880		11.7898
REL STD DEV	1.379		0.342
H2O adjust (mg/l*10k)	461		365

Atmospheric Pressure = 956

*****CALIBRATION SUCCESSFUL*****