

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 4198

A. General Setup and Checks:

1. Initial start up
 - Power On
 - 5-minute countdown during initial warm-up
 - 1-minute countdown to begin test after standby mode
2. Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 - Level 1 (Operator)
 - Level 2 (Field Inspector)
 - Level 3 (Lab Personnel)
4. Printer Setup (Level 3,S,P)
5. Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

1. Internal Printer Test (Level 1,P). Sign and attach record.
 - Paper Feed
2. External Printer Test (Level 1,P). Sign and attach record.
3. Barcode Reader (Level 1,B) (Tone sounds after each)
 - Operator ID
 - Invalid ID
 - ND License
 - R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

1. Continuous Air Blank (Level 2,A)
2. Diagnostic (Level 2,D). Sign and attach test record.
3. Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
Inst Ser#: 80-004198 DSP Ver: 1362.14
Core Ver: 8104.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 12680 s: 1 C: 47.00
9: 12662 s: 2 B: 45

DVM (3 & 9) between 9,000-16,000 counts

RFI indicator with keyed radio: 0 → 1

2. Function (F) Keys:

F5 – Tone

F9 - Fan

F6 – Green LED

F10 – Pump

F7 – Red LED

F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 948 CF: 1.06

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 222 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 47.5 I: 0.0

B: 45.1 V: 0.0

Cell Temperature (C) between 47 ± 0.02 °C

Breath Hose Temperature (B) between 35.5 – 47 °C

Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 155 psi Regulator: 150 ²⁵⁰ psi ^{AEN 1/30/25}

Display and Regulator ± 50 psi of each other Yes or No

Gas tank tare necessary? Yes or No ^{AEN 01/30/25}

If Yes, display readings after tare (Level 3,M,C,G):

Display: 238 psi Regulator: 250 psi

7. General Setup (Level 3,S,G)

8. Configure Start Test (Level 3,S,T)

9. Network (Level 3,S,N)

E. Tests (Sign and attach all test records):

1. Configure wet bath simulator for the following test (Level 1,S)

Zero AC-Wet Calibration Check (Level 1,C):

Sim ser #: MP3063 Lot #: NIA Exp date: NIA

Known value: 0.000 AC

Results are Zero AC

2. Configure wet bath simulator for the following test (Level 1,S)

Low AC-Wet Calibration Check (Level 1,C):

Known Value ≤ 0.03 AC: 0.020 AC

Sim ser #: MP6037 Lot #: 202311D Exp date: 28Nov25

Results ± 0.005 of known AC

3. Configure wet bath simulator for the following test (Level 1,S)
 High AC-Wet Calibration Check (Level 1,C):
 Known Value \geq 0.25 AC: 0.300 AC
 Sim ser #: MP3048 Lot #: 202402C Exp date: 14Feb26
 Results \pm 5% of known AC
4. Configure dry gas standard for the following tests (Level 1,S)
 Known Value: 0.080 AC
 Gas Cylinder Lot #: 01923080A3 Cylinder #: 22 Exp. Date: 2/5/2025
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):
 Sim. Ser #: DR7352 Lot #: ICS 8 Exp. Date: NIA
 Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B):
 Sim. Ser #: DR7111 Lot #: ACE4 Exp. Date: NIA
 Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):
 Sim. Ser #: DR7345 Lot #: MET4 Exp. Date: NIA
 Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):
 Sim. Ser #: DR7350 Lot #: IPA3 Exp. Date: NIA
 Display Reads "Interferent Detect"
9. Ambient Fail (Level 1,B): Display read and test record correct.
10. Invalid Sample (Level 1,B): Display read and test record correct.
11. Sequence Aborted (Level 1,B): Display read test record correct.
12. ABA (Level 1,B): Test record correct.
13. Calibration check out of tolerance (CMS mode): Display read and test record correct.
14. RFI Detected (B,C, or CMS mode): Display read and test record correct.
15. Tank Pressure Below Minimum (CMS mode): Display reads correct.
16. Deficient Samples (CMS mode): Display read and test record correct.
17. Improper Sample (CMS mode): Display read and test record correct.
18. Difference Too Great (CMS mode): Display read and test record correct.
19. CMS (CMS mode): Test record correct.
20. Stability Test (Level 3,T,S):
 Ave: 0.0800 Std Dev: 0.0005 Rel Std Dev: 0.0496

Maintenance/Repairs/Notes: NIA

Instrument is acceptable to be used in the field. Yes or

No

If No, state reason(s) why:

This instrument has had multiple issues with not turning on,
so this instrument will be a classroom only instrument.

Alan Muhl
Inspector Signature

30 January 2025
Date

Janelle Pntschiller
Reviewer

31 Jan. 2025
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 13:10

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNQRSTU

24 characters per line
abcdefghijklmnopqrstuvw
ABCDEFGHIJKLMNQRSTUUVW

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNQRSTUWXYZ123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXYZ1234567890ABCDEF

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On MemFull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 28423080A3
Standard Cyl #? 37
Standard Expiration? 11/05/2025
Oper No? 133237

Flow Cal. Date: 01/09/2012
Slope 637
Intercept -467208

IR Calibration Date: 05/18/2016
3um 9um

0th Coef(*100): -19233 -23050
1st Coef(*100): 267031 132119
2nd Coef(*100): 2449 1142
H2O adj(mg/l*10k): 329 398

***** Test End *****


Operator Signature
ANNA NAREHOOD

Remarks:
Internal Print Test

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 13:28

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890- = |
ABCDEFGHIJKLMNopqrstuvwxyz!@#%\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 28423080A3
Standard Cyl #? 37
Standard Expiration? 11/05/2025
Oper No? 133237

Flow Cal. Date: 01/09/2012
Slope 637
Intercept -467208

IR Calibration Date: 05/18/2016

	3um	9um
0th Coef(*100):	-19233	-23050
1st Coef(*100):	267031	132119
2nd Coef(*100):	2449	1142
H2O adj(mg/l*10k):	329	398

***** Printer Test End *****



Operator Signature
ANNA NAREHOOD

Remarks: External Print Test


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 13:26

Test	AC	Time
01 Room Air	0.000	13:26
02 *Subject Test	REF*	13:27
03 Room Air	0.000	13:27

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = REFUSAL
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 37
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Refusal*

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 13:28

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Diagnostic*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 13:35

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:36
02 Std. Sol.	0.000	13:37
03 Room Air	0.000	13:37
04 Std. Sol.	0.000	13:38
05 Room Air	0.000	13:38
06 Std. Sol.	0.000	13:39
07 Room Air	0.000	13:40

08 Sim Temp = 34.0°C

Simul Ser No = MP3063

Std Sol No = NA

County = 08

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Zero AC calibration check
(0.000 AC)

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 13:41

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:42
02 Std. Sol.	0.020	13:43
03 Room Air	0.000	13:43
04 Std. Sol.	0.020	13:44
05 Room Air	0.000	13:45
06 Std. Sol.	0.021	13:45
07 Room Air	0.000	13:46

08 Sim Temp = 34.0°C

Simul Ser No = MP6037
Std Sol No = 202311D
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Low Ac check*
0.020AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:02

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:03
02 Std. Sol.	0.296	14:04
03 Room Air	0.000	14:04
04 Std. Sol.	0.298	14:05
05 Room Air	0.000	14:06
06 Std. Sol.	0.298	14:06
07 Room Air	0.000	14:07


08 Sim Temp = 34.0°C

Simul Ser No = MP3048

Std Sol No = 202402C

County = 08

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: High AC Check
0.300AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:12

Test	AC	Time
01 Room Air	0.000	14:13
02 *Subject Test	INT*	14:13
03 Room Air	0.000	14:14

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10% Ethanol + 0.05% Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:16

Test	AC	Time
01 Room Air	0.000	14:16
02 *Subject Test	INT*	14:17
03 Room Air	0.000	14:17

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: 0.05% Acetone
Interferent check

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:18

Test	AC	Time
01 Room Air	0.000	14:18
02 *Subject Test	INT*	14:19
03 Room Air	0.000	14:20

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10% Methanol

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:35

Test	AC	Time
01 Room Air	0.000	14:36
02 *Subject Test	INT*	14:37
03 Room Air	0.000	14:37

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10% Isopropanol

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:50

Test	AC	Time
01 Room Air	AMB*	14:50
02 Room Air	0.000	14:51

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = AMBIENT
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Ambient Fail*

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:55

Test	AC	Time
01 Room Air	0.000	14:56
02 Invalid Sample	X.XXX	14:56
03 Room Air	0.000	14:57

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INVALID SAMPLE
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Invalid Sample

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:58

Test	AC	Time
01 Room Air	0.000	14:59
02 *Subject Test	ABT*	14:59
03 Room Air	0.000	14:59

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = SEQUENCE ABORTE
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Sequence Aborted.*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 14:59

Test	AC	Time
01 Room Air	0.000	15:00
02 Subject Test 1	0.000	15:00
03 Room Air	0.000	15:01
04 Reported AC	0.000	15:00

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = ABA TEST
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: ABA test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 15:01

Test	AC	Time
01 Diagnostic	OK	15:02
02 Room Air	0.000	15:02
03 Subject Test 1	0.000	15:03
04 Room Air	0.000	15:05
05 Std. Gas	0.000*	15:06
06 Room Air	0.000	15:07

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = CAL CHECK OOT
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: calibration check
out of tolerance

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 15:07

Test	AC	Time
01 Diagnostic	OK	15:08
02 Room Air	0.000	15:08
03 *Subject Test	RFI*	15:08
04 Room Air	0.000	15:09

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = RFI TEST
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: RFI Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 15:09

Test	AC	Time
01 Diagnostic	OK	15:10
02 Room Air	0.000	15:11
03 *Subject Test	0.000*	15:14
04 Room Air	0.000	15:17
05 Std. Gas	0.079	15:18
06 Room Air	0.000	15:19
07 *Subject Test	0.000*	15:22
08 Room Air	0.000	15:23


No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = DEFICIENT
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Deficient Samples

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 15:24

Test	AC	Time
01 Diagnostic	OK	15:24
02 Room Air	IPS*	15:25
03 Room Air	0.000	15:25

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = IMPROPER SAMPLE
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Improper Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/29/2025 15:49

Test	AC	Time
01 Diagnostic	OK	15:50
02 Room Air	0.000	15:50
03 Subject Test 1	0.000*	15:51
04 Room Air	0.000	15:53
05 Std. Gas	0.080	15:55
06 Room Air	0.000	15:56
07 Subject Test 2	0.393*	15:56
08 Room Air	0.000	15:57

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Difference Too Great

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/30/2025 08:15


Test	AC	Time
01 Diagnostic	OK	08:17
02 Room Air	0.000	08:17
03 Subject Test 1	0.000	08:17
04 Room Air	0.000	08:19
05 Std. Gas	0.079	08:21
06 Room Air	0.000	08:22
07 Subject Test 2	0.000	08:22
08 Room Air	0.000	08:23
09 Reported AC	0.000	08:17

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = CMS TEST
Dr. Lic. = ND/DIS821456
Lot No = 01923080A3
Cyl No = 22
Expiration Date = 02/05/2025
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: CMS Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004198
Location = TOXL 8164.14.00 09/16
01/30/2025 08:24

Test	AC	Time
01 Room Air	0.000	08:25
02 Std. Gas	0.079	08:25
03 Room Air	0.000	08:26
04 Std. Gas	0.079	08:26
05 Room Air	0.000	08:26
06 Std. Gas	0.081	08:27
07 Room Air	0.000	08:27
08 Std. Gas	0.080	08:28
09 Room Air	0.000	08:28
10 Std. Gas	0.080	08:28
11 Room Air	0.000	08:29
12 Std. Gas	0.080	08:29
13 Room Air	0.000	08:30
14 Std. Gas	0.080	08:30
15 Room Air	0.000	08:31
16 Std. Gas	0.080	08:31
17 Room Air	0.000	08:31
18 Std. Gas	0.080	08:32
19 Room Air	0.000	08:32
20 Std. Gas	0.080	08:33
21 Room Air	0.000	08:33
22 Std. Gas	0.080	08:34
23 Room Air	0.000	08:34
24 Std. Gas	0.080	08:34
25 Room Air	0.000	08:35
26 Std. Gas	0.079	08:35
27 Room Air	0.000	08:36
28 Std. Gas	0.080	08:36
29 Room Air	0.000	08:37
30 Std. Gas	0.080	08:37
31 Room Air	0.000	08:37
32 Std. Gas	0.081	08:38
33 Room Air	0.000	08:38
34 Std. Gas	0.080	08:39
35 Room Air	0.000	08:39
36 Std. Gas	0.080	08:39
37 Room Air	0.000	08:40
38 Std. Gas	0.080	08:40
39 Room Air	0.000	08:41
40 Std. Gas	0.080	08:41
41 Room Air	0.000	08:42
42 Cal Check Stats:		
43 Average	0.0800	
44 Std Dev	0.0005	
45 Rel Std Dev	0.64%	

Am Miel

Operator Signature
ANNA NAREHOOD

Remarks: Stability Test

Form 106-I8000