

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

	30			
Reason for Install/Repair				
Install After Receiving From Crime Laboratory Install After Location Change				
Other (Specify)				
Check When Done:				
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
3. Breath Tube Heated.				
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
💢 6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	ring Test).			
7. Repair and/or Maintenance Performed (if any):				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ncy.			
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	heckout (SFN59281, Form			
Field Inspector Signature Alexis Peterson - a. pla 4402	Date 02/11/25			
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to deter the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crim County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.				
Reviewed/Certified By Anna Narehood	Certified Date			

18 Feb 2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006490
Location = ASHL 8164.14.00 09/16
02/11/2025 14:50

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 01524080A4 20 03/05/2026 Standard Expiration? Oper No? 132828

Flow Cal. Date: 09/04/2020

Slope 649 Intercept -598817

IR Calibration Date: 09/04/2020 3um 9um

Oth Coef(*100): -25190 -16392 1st Coef(*100): 271645 136906 2nd Coef(*100): 2533 1222 H2O adj(mg/l*10k): 460 309

******* Printer Test End ********

Operator Signature ALEXIS PETERSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006490 Location = ASHL 8164.14.00 09/16 02/11/2025

14:51

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	14:52
02	Std.	Gas	0.081	14:52
03	Room	Air	0.000	14:53
04	Std.	Gas	0.082	14:53
05	Room	Air	0.000	14:54
	Std.		0.081	14:54
	Room		0.000	14:55

Lot No = 01524080A4

Cyl No = 20

Exp Date = 03/05/2026

County = 26

Oper No. = 132828

Operator Signature ALEXIS PETERSON

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = ASHL 8164.14.00 09/16 02/11/2025

SN 80-006490 14:55

Test	AC	Time
01 Room Air	RFI*	14:56
02 Room Air	0.000	14:57

*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 02/11/2000

Sub Sex = Female Weight = 120 Cit = 0000Test = OTH

Dr. Lic. = ND/0000Lot No = 01524080A4

Cyl No = 20

Expiration Date = 03/05/2026

County = 26

Oper No. = 132828

Operator Signature ALEXIS PETERSON

Remarks:

Form 106-I8000