

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 -00 5 9 45	Instrument Location			
Reason for Install/Repair				
☐ Install After Receiving From Crime Laboratory ☐ Install Afte	r Location Change			
Other (Specify)				
Check When Done:	à			
1. Surge Protector Installed/Property Grounded.				
2. Telephone Line Connected to Intoxilyzer® 8000.				
∑ 3. Breath Tube Heated.				
∠ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
ズ 6. Run Tests:				
A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).				
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C, Key Radio During Test).				
7. Repair and/or Maintenance Performed (if any):				
∑ 8. Complete the Top Portion of the Intoxilyzer® Record (SFN504)	496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.		
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G)	at the Intoxilyzer® Location at the Age	ency.		
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	kilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form		
Field Inspector Signature		Date		
In Blet		3-3-25		
Crime Laboratory Use Only				
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.				
Reviewed/Certified By	An I A	Certified Date		
Anna Narehood (/www.hill.com/ 04March2				

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = BSPD 8164.14.00 09/16
03/03/2025 15:13

********* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 14323080A4 Standard Cyl #? Standard Expiration? 06/05/2025 Oper No? 069823

Flow Cal. Date: 06/20/2013 Slope 637 Intercept -566407

IR Calibration Date: 05/22/2015
3um 9um

0th	Coef(*100):	-21213	-23093
1st	Coef(*100):	279171	139292
2nd	Coef(*100):	2898	1291
H20	adj (mg/l*10k):	332	319

******* Printer Test End *******

Operator Signature TIMOTHY BLETH

Remarks: Print Test

New Install - moved rooms

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BSPD
03/03/2025

Alcohol Analyzer SN 80-005945 8164.14.00 09/16 15:16

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:17
02 Std. Gas	0.081	15:17
03 Room Air	0.000	15:18
04 Std. Gas	0.080	15:18
05 Room Air	0.000	15:18
06 Std. Gas	0.081	15:19
07 Room Air	0.000	15:19

Lot No = 14323080A4

Cyl No = 8

Exp Date = 06/05/2025

County = 08

Oper No. = 069823

Operator Signature TIMOTHY BLETH

Remarks: ACA-QK

New Install - Moved rooms

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BSPD
03/03/2025

Alcohol Analyzer SN 80-005945 8164.14.00 09/16 15:14

DRY CAL CHECK

Test		AC	Time
01 Room	n Air	RFI*	15:14
02 Room	n Air	0.000	15:15

*Invalid Test Inhibited - RFI

Lot No = 14323080A4

Cyl No = 8

Exp Date = 06/05/2025

County = 08

Oper No. = 069823

Operator Signature TIMOTHY BLETH

To Ble

Remarks: RFI Test - OK

New Install - Moved rooms

Form 106-I8000