

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number (A) A A (((()))) Instru	
80-00099/	ment Location alls County Sheriff Office
Reason for Install/Repair	
Install After Receiving From Crime Laboratory Install After Local	tion Change
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Lev	el 2, Function E).
.5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
✓ 6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1	, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any):	,
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Fo	orm 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the	Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzero 104-G), Print Test, ACA Test, and RFI Test.	
Field Inspector Signature Suple Playar	Date 03-10-2020
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be used f the date the Field Inspector performed the installation. This record on file at the County of Burleigh, North Dakota, is certified to be a true and correct copy of the	Office of Attorney General Crime Laboratory Division in the
Reviewed/Certified By Anna Narehood	Certified Date

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006491
Location = WESO 8164.14.00 09/16
03/10/2025 13:53

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 14323080A1 Standard Cyl #? 19 Standard Expiration? 06/05/2025 Oper No? 133682

Flow Cal. Date: 05/22/2018 Slope 639 Intercept -618311

IR Calibration Date: 05/23/2018

3um 9um

Oth Coef(*100): -23641 -17066

1st Coef(*100): 264855 137703

2nd Coef(*100): 2663 1471

H2O adj(mg/l*10k): 556 451

****** Printer Test End *******

Operator Signature
TYLER RISPA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WESO
03/10/2025

Alcohol Analyzer SN 80-006491 8164.14.00 09/16 13:54

DRY CAL CHECK

Test	AC	Time	
01 Room Air	0.000	13:54	
02 Std. Gas	0.080	13:55	
03 Room Air	0.000	13:55	
04 Std. Gas	0.080	13:56	
05 Room Air	0.000	13:56	
06 Std. Gas	0.080	13:57	
07 Room Air	0.000	13:57	

Lot No = 14323080A1

Cyl No = 19

Exp Date = 06/05/2025

County = 52

Oper No. = 133682

Operator Signature TYLER RISPA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = WESO 03/10/2025

SN 80-006491 8164.14.00 09/16 13:58

DRY CAL CHECK

Test			AC	Time
01	Room	Air	RFI*	13:59
02	Room	Air	RFI*	13:59

*Invalid Test Inhibited - RFI

Lot No = 14323080A1

Cyl No = 19

Exp Date = 06/05/2025

County = 52

Oper No. = 133682

Operator Signature TYLER RISPA

Remarks:

Form 106-I8000