

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Seria	al Number 80 -00 4188	Instrument Location	E PD		
	on for Install/Repair				
	Install After Receiving From Crime Laboratory Install After Location Change				
[Other (Specify)				
Chec	k When Done:				
[Surge Protector Installed/Property Grounded.				
Į	2. Telephone Line Connected to Intoxilyzer® 8000.				
I	3. Breath Tube Heated.				
[4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).				
[5. Scan/Enter Gas Cylinder Information (Level 1, Function S).				
[6. Run Tests:				
==,	A. Print Test (Level1, Function P).				
	B. ACA Test (Level 1, Function C).				
ÆΝ	C. Radio Frequency Interference (RFI) Test (CMS Mode or	r Level 1, Function B or C; Key Radio De	uring Test).		
	7. Repair and/or Maintenance Performed (if any):	-			
]	8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.				
[9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	G) at the Intoxilyzer® Location at the Ag	ency.		
Į.	10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.				
Field	Inspector Signature A	HALL 133223	Date 4-4-25		
0:	e Laboratory Use Only				
	•	and the state of t			
the d	installation has been reviewed and the instrument is approved to b ate the Field Inspector performed the installation. This record on f ity of Burleigh, North Dakota, is certified to be a true and correct co	ile at the Office of Attorney General, Cri			
Revie	ewed/Certified By	An I A	Certified Date		
	Anna Narehood (// // // //	Intel	16April2025		

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = KMPD 8164.14.00 09/16
04/04/2025 18:22

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Yes Disable On Memfull? # of Print Copies? Select Std (D/W/I)? Dry 0.08 Standard Value? 0.080 Standard Lot #? 01524080A4 Standard Cyl #? 26 Standard Expiration? 03/05/2026 133223 Oper No?

Flow Cal. Date: 04/30/2020 Slope 645 Intercept -637424

IR Calibration Date: 05/06/2020 3um 9um

Oth Coef(*100): -30092 -33449

1st Coef(*100): 293453 143538

2nd Coef(*100): 151 602

H2O adj(mg/l*10k): 425 462

****** Printer Test End ********

Operator Signature

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = KMPD
04/04/2025

Alcohol Analyzer SN 80-004188 8164.14.00 09/16 18:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:24
02 Std. Gas	0.082	18:24
03 Room Air	0.000	18:24
04 Std. Gas	0.083	18:25
05 Room Air	0.000	18:25
06 Std. Gas	0.082	18:26
07 Room Air	0.000	18:26

Lot No = 01524080A4

Cyl No = 26

Exp Date = 03/05/2026

County = 51

Oper No. = 133223

Operator Signature
ALLEN HALL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = KMPD
04/04/2025

SN 80-004188 8164.14.00 09/16 18:28

DRY CAL CHECK

Test	AC RFI*	Time
01 Room Air		
02 Room Air	RFI*	18:29

*Invalid Test Inhibited - RFI

Lot No = 01524080A4

Cyl No = 26

Exp Date = 03/05/2026

County = 51

Oper No. = 133223

Operator Signature ALLEN HALL

Remarks:

Form 106-I8000