

## **INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Instrument Location		
Reason for Install/Repair	* e	
│		
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
1 2. Telephone Line Connected to Intoxilyzer® 8000.		
⅓ 3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).		
✓ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
∑. 6. Run Tests:		
🔀 A. Print Test (Level1, Function P).		
☑ B. ACA Test (Level 1, Function C).		
☑ C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio D	Ouring Test).	
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the	Intoxilyzer® for Use.	
120-G) at the Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the A	gency.	
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair 104-G), Print Test, ACA Test, and RFI Test.	Checkout (SFN59281, Form	
Field Inspector Signature	Date ,	
1///	4/7/2025	
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to det the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.		
Reviewed/Certified By	Certified Date	
Anna Narehood	16 A pril 2025	

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007097
Location = BURK 8164.16.00 09/18
04/07/2025 11:05

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Yes Display Third Digit? Inhib Printer(Y/N)? No Display Volume? No Yes Disable On Memfull? # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Value? Standard Lot #? 25024080A3 Standard Cyl #? 57 Standard Expiration? 10/05/2026 Oper No? 132664

Flow Cal. Date: 09/03/2020 Slope 660 Intercept -636264

IR Calibration Date: 04/04/2024
3um 9um

Oth Coef(\*100): -25634 -25706
1st Coef(\*100): 238131 134910
2nd Coef(\*100): 3596 1460
H2O adj(mg/l\*10k): 614 521

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature CHRISTOPHER ALMLIE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BURK
04/07/2025

Alcohol Analyzer SN 80-007097 8164.16.00 09/18 11:06

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:07
02 Std. Gas	0.079	11:08
03 Room Air	0.000	11:08
04 Std. Gas	0.078	11:09
05 Room Air	0.000	11:09
06 Std. Gas	0.079	11:09
07 Room Air	0.000	11:10

Lot No = 25024080A3

Cyl No = 57

Exp Date = 10/05/2026

County = 07

Oper No. = 132664

Operator Signature CHRISTOPHER ALMLIE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007097 Location = BURK 8164.16.00 09/18 04/07/2025

12:02

Test			AC	Time
01	Room	Air	RFI*	12:03
02	Room	Air	0.000	12:04

\*Invalid Test Inhibited - RFI

Sub Name = TEST, TEST TEST

Sub DOB = 04/07/1991

Sub Sex = Male

Weight = 100Cit = NA

Test = DUI

Dr. Lic. = ND/NA

Lot No = 25024080A3

Cyl No = 57

Expiration Date = 10/05/2026

County = 07

Oper No. = 132664

Operator Signature CHRISTOPHER ALMLIE

Remarks:

Form 106-I8000