

SFN 59281 (06/2018)

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

Serial Number 00 42-00	Instrument Location		
Reason for Install/Repair			
🖄 Install After Receiving From Crime Laboratory 🛛 🗌 Install Afte	er Location Change		
Other (Specify)			
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
X_2. Telephone Line Connected to Intoxilyzer® 8000.			
🖾 3. Breath Tube Heated.			
🗹 4. Enter Preliminary Data (I.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).		
🕅 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
A 6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).	
🔀 7. Repair and/or Maintenance Performed (if any):	14L DSP		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature		Date 4-17-2025	

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By Certified Date Anna Narehood 23April2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004200 Location = NAPO 8164.14.00 09/16 04/17/2025 14:46 *********** Printer Test *********** abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Prelim KSIC. Display Third Digit? Yes Display Third Digit? No Yes Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? 25024080A3 Standard Cyl #? 60 Standard Expiration? 10/05/2026 131100 Oper No? Flow Cal. Date: 06/25/2019 614 Slope -551153 Intercept IR Calibration Date: 06/25/2019 3um 9um ______ 0th Coef(*100):-187921st Coef(*100):2670012nd Coef(*100):1903 -22164 140777 1083 H2O adj(mg/l*10k): 403 445 *********** Printer Test End **********

M

Operator Signature MICHAEL MAROLT

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004200
Location = NAPO	8164.14.00 09/16
04/17/2025	14:47

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.082 0.000	14:48 14:48 14:49 14:49 14:49 14:50

Mil le

Operator Signature MICHAEL MAROLT

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer SN 80-004200 North Dakota Model 8000 8164.14.00 09/16 Location = NAPO04/17/2025 14:51 Test Time AC 14:53 01 Room Air RFI* 14:53 02 Room Air 0.000 *Invalid Test Inhibited - RFI Sub Name = TEST, TEST TEST Sub DOB = 12/20/1994Sub Sex = Male Weight = 500Test = DUI Cit = 11111111111Dr. Lic. = ND/TES946541Lot No = 25024080A3Cyl No = 60Expiration Date = 10/05/2026Oper No. = 131100County = 24

MA MIL

Operator Signature MICHAEL MAROLT

Remarks:

Form 106-18000