

## **INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59281 (06/2018)

Serial Number SG - COG 494	Instrument Location	
Reason for Install/Repair		
X Install After Receiving From Crime Laboratory	Install After Location Change	
Other (Specify)		

Check When Done:

K)	1.	Surge	Protector	Installed/Pr	roperty	Grounded.
----	----	-------	-----------	--------------	---------	-----------

X 2. Telephone Line Connected to Intoxilyzer® 8000.

X 3. Breath Tube Heated.

X 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).

5. Scan/Enter Gas Cylinder Information (Level 1, Function S).

A 6. Run Tests:

A. Print Test (Level1, Function P).

B. ACA Test (Level 1, Function C).

🕅 C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).

7. Repair and/or Maintenance Performed (if any): M/A
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

1 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.

10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature	Date
1 Ada	4-15-2025

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

	$ \sim$	)	A			$\leq$	
Reviewed/Certified By	A	/ /		$\Lambda \parallel$		/	Certified Date
Anna Narehood	Λ	m	/	Inh	X		29April2025
	_						

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006494 Location = CAVA 8164.14.00 09/16 23:41 04/19/2025 \*\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*\*\*\*\* abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+? Current Instrument Setup Enabled Data Entry Mode: Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? Yes Disable On Memfull? # of Print Copies? 1 Dry 0.08 Select Std (D/W/I)? 0.080 Standard Value? 25024080A3 Standard Lot #? Standard Cyl #? 007 Standard Expiration? 10/05/2026 020832 Oper No? 08/18/2015 Flow Cal. Date: 655 Slope -624238 Intercept IR Calibration Date: 05/10/2017 9um 3um \_\_\_\_\_ Oth Coef(\*100): -28969 1st Coef(\*100): 256589 2nd Coef(\*100): 2527 -21478 140860 1330 473 H2O adj(mg/l\*10k): 655 \*\*\*\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*\*\*

Operator Signature CORY MORTENSEN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006494Location = CAVA8164.14.00 09/1604/19/202523:45

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.081 0.000 0.081 0.000 0.081 0.000	23:45 23:46 23:46 23:46 23:47 23:47 23:47

Lot No = 25024080A3 Cyl No = 007 Exp Date = 10/05/2026 County = 34 Oper No. = 020832

#4505 Operator Signature

CORY MORTENSEN

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006494Location = CAVA8164.14.00 09/1604/19/202523:49

DRY	CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Room Air	0.000 0.082 RFI* RFI*	23:49 23:50 23:50 23:50 23:50
*Invalid Test Inhibited - RH	ľ	

Operator Signature CORY MORTENSEN

Remarks:

Test 27

Form 106-18000