

## **INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-064945	Instrument Location GFTN				
Reason for Install/Repair					
Install After Receiving From Crime Laboratory ☐ Install After Location Change					
Other (Specify)					
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
2. Telephone Line Connected to Intoxilyzer® 8000.					
3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
AEN 🔀 6. Run Tests:					
A. Print Test (Level1, Function P).					
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).			
7. Repair and/or Maintenance Performed (if any):					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0					
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.					
Field Inspector Signature		Date 4-78-2025			
Crime Laboratory Use Only					
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.					
Reviewed/Certified By Anna Narehood	M 10	Certified Date			
Anna Narehood //www / Land 29April2025					

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004945
Location = GFTN 8164.14.00 09/16
04/28/2025 07:30

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01524080A4 Standard Cyl #? 31 Standard Expiration? 03/05/2026 Oper No? 133368

Flow Cal. Date: 09/26/2012 Slope 622 Intercept -485836

IR Calibration Date: 01/21/2025
3um 9um

Oth Coef(\*100): -26009 -27015
1st Coef(\*100): 264294 142861
2nd Coef(\*100): 2536 1231
H2O adj(mg/l\*10k): 461 365

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature
JOHN GROSSBAUER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFTN
04/28/2025

Alcohol Analyzer SN 80-004945 8164.14.00 09/16 07:31

## DRY CAL CHECK

Те	est		AC	Time
01	Room	Air	0.000	07:31
02	Std.	Gas	0.078	07:32
03	Room	Air	0.000	07:32
04	Std.	Gas	0.078	07:33
05	Room	Air	0.000	07:33
06	Std.	Gas	0.079	07:33
07	Room	Air	0.000	07:34

Lot No = 01524080A4

Cyl No = 31

Exp Date = 03/05/2026

County = 50

Oper No. = 133368

Operator Signature JOHN GROSSBAUER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFTN
04/28/2025

Alcohol Analyzer SN 80-004945 8164.14.00 09/16 07:34

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	RFI*	07:35
02	Room	Air	RFI*	07:35

\*Invalid Test Inhibited - RFI

Lot No = 01524080A4

Cyl No = 31

Exp Date = 03/05/2026

County = 50

Oper No. = 133368

Operator Signature JOHN GROSSBAUER

Remarks:

Form 106-I8000