

## **INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-00667/	Instrument Location DEVL
Reason for Install/Repair	Neve
Install After Receiving From Crime Laboratory	nstall After Location Change
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), an	nd Location; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Fund	ction S).
☑-6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS)	Mode or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any):	
	d (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
△9. File Previous Intoxilyzer® Record (SFN504096, Formula 1)  1. Second (SFN504096, Formula 2)  2. Second (SFN504096, Formula 2)  2. Second (SFN504096, Formula 2)  3. Second (SFN504096, Formula 2)  3. Second (SFN504096, Formula 2)  4. Seco	rm 120-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Complete 104-G), Print Test, ACA Test, and RFI Test.	oleted Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
Field Inspector Signature	Date
	14534 04/28/202
Crime Laboratory Use Only	
	oved to be used for the analysis of breath to determine alcohol concentration fro cord on file at the Office of Attorney General, Crime Laboratory Division, in the correct copy of the documents received.
Reviewed/Certified By	Certified Date
Anna Narehood /	// 30April2025

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006671 Location = DEVL 8164.14.00 09/16 04/28/2025 20:51

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefqhijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 25024080A3 Standard Cyl #? 4 Standard Expiration? 10/05/2026 Oper No? 132008

Flow Cal. Date: 08/22/2016 Slope 703 Intercept -691868

IR Calibration Date: 08/22/2016 3um 9um

Oth Coef(\*100): -17135 1st Coef(\*100): 254798 2nd Coef(\*100): 3386 -17135 -24345 132856 1354 H20 adj(mq/l\*10k): 610

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature

WILLIAM EARL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DEVL
04/28/2025

Alcohol Analyzer SN 80-006671 8164.14.00 09/16 20:58

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	20:59
02	Std.	Gas	0.082	20:59
03	Room	Air	0.000	21:00
04	Std.	Gas	0.082	21:00
05	Room	Air	0.000	21:01
06	Std.	Gas	0.082	21:01
07	Room	Air	0.000	21:01

Lot No = 25024080A3

Cyl No = 4

Exp Date = 10/05/2026

County = 36

Oper No. = 132008

Operator Signature WILLIAM EARL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006671 Location = DEVL 8164.14.00 09/16 04/28/2025

21:03

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:04
02 Std. Gas	0.083	21:05
03 Room Air	RFI*	21:05
04 Room Air	0.000	21:05

\*Invalid Test Inhibited - RFI

Lot No = 25024080A3

Cyl No = 4

Exp Date = 10/05/2026

County = 36

Oper No. = 132008

Operator Signature WILLIAM EARL

Remarks:

Form 106-I8000