

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location		
Reason for Install/Repair	Killdeer P.D.		
Install After Receiving From Crime Laboratory Install After	er Location Change		
Other (Specify)	Ŭ		
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
☑ 3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	ion; Level 2, Function E).		
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	ring Test).	
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) and Place it by the li	ntoxilyzer® for Use.	
☑ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	6) at the Intoxilyzer® Location at the Age	ncy.	
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature		Date	
P.MM		4/24/25	
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.			
Reviewed/Certified By Anna Narehood	Mall	Certified Date 07May2025	

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = KDPD

04/29/2025

Alcohol Analyzer SN 80-004938 8164.14.00 09/16 14:54

****** Printer Test *******

12 char/line abcdefghijkl ABCDEFGHIJKL

16 char per line abcdefghijklmnop ABCDEFGHIJKLMNOP

21 chars per line abcdefghijklmnopqrstu ABCDEFGHIJKLMNOPQRSTU

24 characters per line abcdefghijklmnopqrstuywx ABCDEFGHIJKLMNOPQRSTUUWX

32 characters per line abcdefghijklmnopqrstuvwxyz123456 ABCDEFGHIJKLMNOPQRSTUUWXYZ123456

42 characters per line abcdefghijklmnopqrstuuwxyz1234567890abcdef ABCDEFGHIJKLMNOPORSTUUWXYZ1234567890ABCDEF

urrent Instrument Se lata Entry Mode: Start Test Sequence: Display Prelim Rslt? Display Third Digit? Unhib Printer(Y/N)? Display Volume? Disable On Memfull? # of Print Copies? Select Std (D/W/1)? Standard Value? Standard Lot #? Standard Cyl #? Standard Expiration Oper No? Flow Cal. Date: Slope Intercept	?	587 -338327	26 20
IR Calibration Date	e: 3um	06/10/20	120 9um
Oth Coef(*100):	-224	34	-2317
lst Coef(*100):	2698	349	13679
2nd Coef(*100):	963		840

******* Test End ********

370

H2O adj(mg/l*10k): 355

Operator Signature PHIL WOLF

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000 SN 80-004938
Location = KDPD 8164.14.00 09/16
04/29/2025 14:43

DRY CAL CHECK

Test ·	AC	Time
01 Room Air	0.000	14:43
02 Std. Gas	0.080	14:44
03 Room Air	0.000	14:44
04 Std. Gas	0.080	14:44
05 Room Air	0.000	14:45
06 Std. Gas	0.079	14:45
07 Room Air	0.000	14:46

Lot No = 25024080A3

Cyl No = 2

Exp Date = 10/05/2026

County = 13

Oper No. = 120547

Operator Signature PHIL WOLF

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = KDPD

04/29/2025

Alcohol Analyzer SN 80-004938 8164.14.00 09/16 14:39

DRY CAL CHECK

Te	25 t	۴	AC	Time
	Room Room	 ***	:[* .000	14:39 14:40

*Invalid Test Inhibited - RFI

Lot No = 25024080A3

Cyl No = 2

Exp Date = 10/05/2026

County = 13

Oper No. = 120547

Operator Signature PHIL WOLF

Remarks:

Form 106-I8000