

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

	1	
Serial Number 80 - 0049 44	Instrument Location OFM	
Reason for Install/Repair		
Install After Receiving From Crime Laboratory 🔲 Install After	er Location Change	
Other (Specify)		
		,
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locat	ion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
☐←C. Radio Frequency Interference (RFI) Test (CMS Mode of	r Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):	A	
28. Complete the Top Portion of the Intoxilyzer® Record (SFN5		ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form
Field Inspector Signature		Date 5/7/2025
I William Com		5/7/2025
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct county.	file at the Office of Attorney General, Crir	
Reviewed/Certified By	M 10	Certified Date
Δnna Narehood / // 🛕	// a / //	-0.014 $-0.00$ $E$

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004944
Location = DEVL 8164.14.00 09/16
05/07/2025 13:19

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 25024080A3 Standard Cyl #? 004 Standard Expiration? 10/05/2026 Oper No? 132008

Flow Cal. Date: 07/01/2011 Slope 685 Intercept -634740

IR Calibration Date: 07/01/2011

Oth Coef(\*100): -22591 -17234

1st Coef(\*100): 265915 132691

2nd Coef(\*100): 2841 1562

H2O adj(mg/l\*10k): 559 425

3um

9um

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature
WILLIAM EARL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = DEVL 05/07/2025

SN 80-004944 8164.14.00 09/16 13:22

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:23
02 Std. Gas	0.082	13:23
03 Room Air	0.000	13:23
04 Std. Gas	0.082	13:24
05 Room Air	0.000	13:24
06 Std. Gas	0.082	13:25
07 Room Air	0.000	13:25

Lot No = 25024080A3

Cyl No = 004

Exp Date = 10/05/2026

County = 36

Oper No. = 132008

Operator Signature WILLIAM EARL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = DEVL 05/07/2025

SN 80-004944 8164.14.00 09/16 13:25

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:26
02 Std. Gas	0.082	13:26
03 Room Air	RFI*	13:26
04 Room Air	0.000	13:27

\*Invalid Test Inhibited - RFI

Lot No = 25024080A3

Cyl No = 004

Exp Date = 10/05/2026

County = 36

Oper No. = 132008

Operator Signature WILLIAM EARL

Remarks:

Form 106-I8000