




INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT
NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59281 (06/2018)

| | |
|---|------------------------------------|
| Serial Number 80-004939 | Instrument Location BCCC |
| Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____ | |


Check When Done:

- ☒ 1. Surge Protector Installed/Property Grounded.
- ☒ 2. Telephone Line Connected to Intoxilyzer® 8000.
- ☒ 3. Breath Tube Heated.
- ☒ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- ☒ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 6. Run Tests:
 - ☒ A. Print Test (Level 1, Function P).
 - ☒ B. ACA Test (Level 1, Function C).
 - ☒ C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- ☐ 7. Repair and/or Maintenance Performed (if any): _____
- ☒ 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

| | |
|--|--------------------------|
| Field Inspector Signature  | Date 5/16/2025 |
|--|--------------------------|

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

| | |
|--|------------------------------------|
| Reviewed/Certified By Anna Narehood  | Certified Date 19May2025 |
|--|------------------------------------|

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = BCCC 8164.14.00 09/16
05/16/2025 15:12

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:13 |
| 02 Std. Gas | 0.080 | 15:13 |
| 03 Room Air | 0.000 | 15:14 |
| 04 Std. Gas | 0.080 | 15:14 |
| 05 Room Air | 0.000 | 15:14 |
| 06 Std. Gas | 0.080 | 15:15 |
| 07 Room Air | 0.000 | 15:15 |

Lot No = 08425080A1
Cyl No = 24
Exp Date = 05/05/2027
County = 02 Oper No. = 100012



Operator Signature
DANA RUSTEBAKKE

Remarks:

ACA Test - Good

Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004939
Location = BCCC 8164.14.00 09/16
05/16/2025 15:16

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 15:16 |
| 02 Std. Gas | RFI* | 15:17 |
| 03 Room Air | 0.000 | 15:17 |

*Invalid Test
Inhibited - RFI

Lot No = 08425080A1
Cyl No = 24
Exp Date = 05/05/2027
County = 02 Oper No. = 100012



Operator Signature
DANA RUSTEBAKKE

Remarks:

RFI Test - Good

Form 106-I8000