

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	
80-005952 Reason for Install/Repair	Wishek PD	
	er Location Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Local	tion; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
☑ B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode of	r Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN	50496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-	G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed In 104-G), Print Test, ACA Test, and RFI Test.	toxilyzer® 8000 Installation and Repair C	Checkout (SFN59281, Form
Field Inspector Signature		Date 5/16/25
Office Laboratory Has Only		
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to the date the Field Inspector performed the installation. This record on County of Burleigh, North Dakota, is certified to be a true and correct of	file at the Office of Attorney General, Crit	
Reviewed/Certified By Anna Narehood	M 10	Certified Date
Allia Naielloud /// /	1 / / /	19May2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005952
Location = WISH 8164.14.00 09/16
05/16/2025 18:23

\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup
Data Entry Mode:

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
# of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 01524080A4
Standard Cyl #? 19
Standard Expiration? 03/05/2026
Oper No? 133473

Flow Cal. Date: 04/16/2021 Slope 657 Intercept -744615

IR Calibration Date: 04/16/2021 3um 9um

Oth Coef(\*100): -24450 -14639
1st Coef(\*100): 266154 132500
2nd Coef(\*100): 2284 1091
H2O adj(mq/l\*10k): 516 323

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature TABITHA KINZLER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005952 Location = WISH 05/16/2025

8164.14.00 09/16 18:24

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:24
02 Std. Gas	0.080	18:25
03 Room Air	0.000	18:25
04 Std. Gas	0.082	18:25
05 Room Air	0.000	18:26
06 Std. Gas	0.081	18:26
07 Room Air	0.000	18:27

Lot No = 01524080A4

Cyl No = 19

Exp Date = 03/05/2026

County = 26

Oper No. = 133473

Operator Signature TABITHA KINZLER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = WISH 05/16/2025

SN 80-005952 8164.14.00 09/16 18:28

## DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	18:28
02 Room Air	0.000	18:29

\*Invalid Test Inhibited - RFI

Lot No = 01524080A4

Cyl No = 19

Exp Date = 03/05/2026

County = 26

Oper No. = 133473

Operator Signature TABITHA KINZLER

Remarks:

Form 106-I8000