

**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59281 (06/2018)

Serial Number 80-005357	Instrument Location GRAF
Reason for Install/Repair <input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

Check When Done:

- ☒ 1. Surge Protector Installed/Property Grounded.
- ☒ 2. Telephone Line Connected to Intoxilyzer® 8000.
- ☒ 3. Breath Tube Heated.
- ☒ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- ☒ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- ☒ 6. Run Tests:
 - ☒ A. Print Test (Level 1, Function P).
 - ☒ B. ACA Test (Level 1, Function C).
 - ☒ C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- ☐ 7. Repair and/or Maintenance Performed (if any): _____
- ☒ 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature <i>Jarvis Miller</i>	Date 5/29/2025
---	--------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By Anna Narehood <i>Anna Narehood</i>	Certified Date 21May2025
---	------------------------------------

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005357
Location = GRAF 8164.14.00 09/16
05/20/2025 17:11

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\$%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 01524080A4
Standard Cyl #? 30
Standard Expiration? 03/05/2026
Oper No? 132597

Flow Cal. Date: 05/16/2018
Slope 573
Intercept 140460

IR Calibration Date: 06/21/2023
3um 9um

0th Coef(*100):	-6039	-13097
1st Coef(*100):	257461	136071
2nd Coef(*100):	3187	1488
H2O adj(mg/l*10k):	315	418

***** Printer Test End *****



Operator Signature
JACOB NIKKO

Remarks:

The Undersigned, having custody of the original public record, certifies that the information hereon is a true and correct copy of the original document maintained as a part of the activities of this agency.

Form 106-I8000

Signature/Title Jacob Nikko/Deputy

Agency/Date Walsh Co. 5/20/25

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005357
Location = GRAF 8164.14.00 09/16
05/20/2025 17:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:21
02 Std. Gas	0.082	17:21
03 Room Air	0.000	17:22
04 Std. Gas	0.081	17:22
05 Room Air	0.000	17:23
06 Std. Gas	0.081	17:23
07 Room Air	0.000	17:24

Lot No = 01524080A4

Cyl No = 30

Exp Date = 03/05/2026

County = 50

Oper No. = 132597



Operator Signature
JACOB NIKKO

Remarks:

Form 106-I8000

The Undersigned, having custody of the
original public record, certifies that the
information hereon is a true and correct
copy of the original document maintained
as a part of the activities of this agency.

Signature/Title Jacob Nikko / Deputy

Agency/Date Walsh Co 50/5/20/2025

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005357
Location = GRAF 8164.14.00 09/16
05/20/2025 17:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:18
02 Std. Gas	0.082	17:18
03 Room Air	0.000	17:19
04 Std. Gas	0.081	17:19
05 Room Air	RFI*	17:19
06 Room Air	0.000	17:20

*Invalid Test
Inhibited - RFI

Lot No = 01524080A4
Cyl No = 30
Exp Date = 03/05/2026
County = 50 Oper No. = 132597



Operator Signature
JACOB NIKKO

Remarks:

Form 106-I8000

The Undersigned, having custody of the
original public record, certifies that the
information hereon is a true and correct
copy of the original document maintained
as a part of the activities of this agency.

Signature/Title Jacob Nikko / Deputy

Agency/Date Walsh CO SO / 5-20-2025