

## **INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number (/ ) / / / / / / / / / / / / / / / / /	
Serial Number 80-106491 Instrument Location WC115 (5) SO	
Reason for Install/Repair	
Install After Receiving From Crime Laboratory 🔲 Install After Location Change	
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
€ 6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	ring Test).
7. Repair and/or Maintenance Performed (if any):	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer®	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	ency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	Checkout (SFN59281, Form
Field Inspector Signature Lylu Market	Date 05/20/202
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to dete the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crit County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	
Reviewed/Certified By Anna Narehood	Certified Date 23May2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006491
Location = WESO 8164.14.00 09/16
05/22/2025 09:23

\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 08425080A1 12 05/05/2027 Standard Expiration? Oper No? 133682

Flow Cal. Date: 05/22/2018 Slope 639 Intercept -618311

IR Calibration Date: 05/23/2018
3um 9um

Oth Coef(\*100): -23641 -17066
1st Coef(\*100): 264855 137703
2nd Coef(\*100): 2663 1471
H2O adj(mg/l\*10k): 556 451

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature TYLER RISPA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WESO
05/22/2025

Alcohol Analyzer SN 80-006491 8164.14.00 09/16 09:23

## DRY CAL CHECK

Test		AC	Time
01 Room	Air	0.000	09:24
02 Std.	Gas	0.080	09:24
03 Room	ı Air	0.000	09:25
04 Std.	Gas	0.079	09:25
05 Room	ı Air	0.000	09:26
06 Std.	Gas	0.079	09:26
07 Room	ı Air	0.000	09:27

Lot No = 08425080A1

Cyl No = 12

Exp Date = 05/05/2027

County = 52

Oper No. = 133682

Operator Signature TYLER RISPA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WESO
05/22/2025

Alcohol Analyzer SN 80-006491 8164.14.00 09/16 09:28

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	RFI*	09:29
02	Room	Air	RFI*	09:29

\*Invalid Test Inhibited - RFI

Lot No = 08425080A1

Cyl No = 12

Exp Date = 05/05/2027

County = 52

Oper No. = 133682

Operator Signature TYLER RISPA

Remarks:

Form 106-I8000