

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

	I
Serial Number	SARCHENT COUNTY SITES IFF WHICE
Reason for Install/Repair	J. J. W. J.
Install After Receiving From Crime Laboratory   Install After	er Location Change
Other (Specify)	
Check When Done:	
Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
A. Print Test (Level1, Function P).  B. ACA Test (Level 1, Function C).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or	
7. Repair and/or Maintenance Performed (if any):	4 Inspection
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	s) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
Field Inspector Signature	Date,
Jon was	Date 5/27/25
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fi County of Burleigh, North Dakota, is certified to be a true and correct co	
Reviewed/Certified By Anna Narehood	Certified Date 30May2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006690
Location = SACO	8164.14.00 09/16
05/27/2025	18:21

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current Instrument Setup	
Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	08425080A1
Standard Cyl #?	11
Standard Expiration?	05/05/2027

Flow Cal. Date:	08/23/2018
Slope	677
Intercept	-593298

IR Calibration Dat	ce:	08/23/2018
	3um	9um
0th Coef(*100):	-8380	-21663
1st Coef(*100):	265142	130429
2nd Coef(*100):	3933	1429
H2O adj(mg/l*10k):	412	539

132607

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature
JONATHAN WOLD

Remarks:

Oper No?

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = SACO
05/27/2025

Alcohol Analyzer SN 80-006690 8164.14.00 09/16 18:29

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:30
02 Std. Gas	0.080	18:30
03 Room Air	0.000	18:31
04 Std. Gas	0.080	18:31
05 Room Air	0.000	18:32
06 Std. Gas	0.080	18:32
07 Room Air	0.000	18:33

Lot No = 08425080A1

Cyl No = 11

Exp Date = 05/05/2027

County = 41

Oper No. = 132607

Operator Signature
JONATHAN WOLD

Jon com

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = SACO
05/27/2025

Alcohol Analyzer SN 80-006690 8164.14.00 09/16 18:33

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:34
02 Reference	RFI*	18:34
03 Room Air	0.000	18:35

\*Invalid Test Inhibited - RFI

Lot No = 08425080A1

Cyl No = 11

Exp Date = 05/05/2027

County = 41

Oper No. = 132607

Operator Signature JONATHAN WOLD

Remarks:

Form 106-I8000