



**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION**

INTOXILYZER® 8000 CALIBRATION ADJUSTMENT

Intoxilyzer® 8000 Serial Number: 80-00 5362 Calibration Adjustment Location: TOXL

A. Pre-Adjustment

Replaced Simulator Return O-Ring Yes or No

B. Calibration Adjustment (Level 3,M,C,O)

1. ☒ Autocalibration Printout Attached
 - ☒ Max Power Res Value ≥ 10
 - ☒ Auto Range Res Value ≥ 4
2. Simulator Solutions for Calibration Adjustment

| Soln. | g/210 L | Lot No. | Exp. Date | Simulator SN |
|-------|---------|-----------------------------|-----------------------------|--------------|
| 1 | 0.000 | NA-Milli-Q H ₂ O | NA-Milli-Q H ₂ O | MP3003 |
| 2 | 0.040 | 202410D | 22Oct2026 | MP6038 |
| 3 | 0.080 | 202501A | 15Jan 2027 | MP3057 |
| 4 | 0.100 | 202408F | 28Sep2026 | MP5319 |
| 5 | 0.300 | 202402C | 14Feb2026 | MP6035 |

3. 0.080 AC Calibration Gas for H₂O Adjustment

Lot No. 14323080A4 Cyl No. 41 Exp. Date: 10/5/25

4. Atmospheric Pressure

Displayed by Intoxilyzer® 8000 935 mbar
Adjusted to using barometer 960 mbar
Auto Calibration Report printout 960 mbar
Barometer Model 10510-922
Barometer Serial Number 250063738
Barometer Calibration Expiration Date 04Feb2027

5. ☒ Screen displayed "Calibration Success"

6. ☒ Calibration Adjustment Printout Attached

- ☒ Solution 1 Avg % Abs ≤ 0.2500
- ☒ Solution 2-5 REL STD DEV ≤ 3.000
- ☒ Residual (g/210 L) values for solutions 1 - 5 ≤ 0.0020 for 3 μ m and 9 μ m channels

☐ Dry Gas H₂O adjustment sum for 3 µm and 9 µm channels within ± 10

3 µm 3393 (Ave.) + 416 (H₂O Adj.) = 3809


9 µm 3371 (Ave.) + 438 (H₂O Adj.) = 3809

C. Is an Annual Inspection due for this instrument? Yes or No

If Yes, complete Intoxilyzer 8000 Annual Inspection (Document ID: 11698)

If No, complete Intoxilyzer 8000 Calibration (Document ID: 11871).

Remarks/Notes: This calibration adjustment was performed due to
high dry gas calibration values during the Annual inspection.



Breath Alcohol Analyst Signature

05 June 2025

Date



Reviewer Signature

05 June 2025

Date

TOXL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005362
06/05/2025 08:37:01

Auto Calibration
Max Power Res Value = 46
Auto Range Res Value = 33

Dr. Muel
05 June 2025

TOXL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005362
06/05/2025 08:37:01

Auto Calibration

pg 1 of 2

| <<<<< 3um >>>>> | | | <<<<< 9um >>>>> | | |
|--|--------|-------------|-----------------|-------------|--|
| ----- | | | | | |
| Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1 | | | | | |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | |
| Sample #1 | 0.0760 | (-0.0210) | 0.1590 | (-0.0220) | |
| Sample #2 | 0.0900 | (0.0180) | 0.1760 | (0.0040) | |
| Sample #3 | 0.0960 | (0.0380) | 0.1750 | (0.0130) | |
| Sample #4 | 0.0350 | (0.0680) | 0.1400 | (0.0220) | |
| Avg % Abs | 0.0737 | (0.0413) | 0.1637 | (0.0130) | |
| STD DEV | 0.0336 | (0.0252) | 0.0205 | (0.0090) | |
| REL STD DEV | 45.639 | (60.886) | 12.527 | (69.231) | |
| ----- | | | | | |
| Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1 | | | | | |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | |
| Sample #1 | 0.8070 | (0.0260) | 1.4940 | (0.0060) | |
| Sample #2 | 0.7800 | (0.0470) | 1.4940 | (0.0190) | |
| Sample #3 | 0.7890 | (0.0520) | 1.4610 | (0.0370) | |
| Sample #4 | 0.8030 | (0.0380) | 1.4710 | (0.0190) | |
| Avg % Abs | 0.7907 | (0.0457) | 1.4753 | (0.0250) | |
| STD DEV | 0.0116 | (0.0071) | 0.0169 | (0.0104) | |
| REL STD DEV | 1.466 | (15.536) | 1.147 | (41.569) | |
| ----- | | | | | |
| Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1 | | | | | |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | |
| Sample #1 | 1.5370 | (-0.0370) | 2.7650 | (-0.0100) | |
| Sample #2 | 1.4890 | (-0.0050) | 2.8060 | (-0.0150) | |
| Sample #3 | 1.4880 | (-0.0120) | 2.7620 | (-0.0250) | |
| Sample #4 | 1.5210 | (-0.0050) | 2.7790 | (-0.0060) | |
| Avg % Abs | 1.4993 | (-0.0073) | 2.7823 | (-0.0153) | |
| STD DEV | 0.0188 | (0.0040) | 0.0222 | (0.0095) | |
| REL STD DEV | 1.252 | (55.111) | 0.797 | (61.985) | |
| ----- | | | | | |
| Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1 | | | | | |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | |
| Sample #1 | 1.9210 | (-0.0300) | 3.4650 | (0.0130) | |
| Sample #2 | 1.8960 | (0.0010) | 3.4880 | (0.0210) | |
| Sample #3 | 1.9010 | (-0.0020) | 3.4670 | (0.0390) | |
| Sample #4 | 1.8440 | (0.0230) | 3.4370 | (0.0410) | |
| Avg % Abs | 1.8803 | (0.0073) | 3.4640 | (0.0337) | |
| STD DEV | 0.0316 | (0.0137) | 0.0256 | (0.0110) | |
| REL STD DEV | 1.679 | (186.142) | 0.740 | (32.718) | |
| ----- | | | | | |
| Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1 | | | | | |
| Sample | % Abs | (% Abs Ref) | % Abs | (% Abs Ref) | |
| Sample #1 | 5.2900 | (-0.0080) | 9.4700 | (-0.0060) | |
| Sample #2 | 5.2280 | (0.0520) | 9.3610 | (0.0900) | |
| Sample #3 | 5.2260 | (0.0620) | 9.3570 | (0.0990) | |
| Sample #4 | 5.2350 | (0.0530) | 9.3770 | (0.0840) | |
| Avg % Abs | 5.2297 | (0.0557) | 9.3650 | (0.0910) | |
| STD DEV | 0.0047 | (0.0055) | 0.0106 | (0.0075) | |
| REL STD DEV | 0.090 | (9.894) | 0.113 | (8.297) | |
| ----- | | | | | |

Handwritten signature
June 2025
AEB

TOXL
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005362
 06/05/2025 08:37:01

Auto Calibration

pg 2 of 2

| <<<< 3um >>>> | | | <<<< 9um >>>> | | |
|-------------------|----------|----------|---------------|----------|----------|
| ----- | | | ----- | | |
| Zero Order Coef | -175.86 | | | -208.11 | |
| First Order Coef | 2581.73 | | | 1386.66 | |
| Second Order Coef | 35.03 | | | 17.16 | |
| ----- | | | ----- | | |
| Act | Fit | Residual | Act | Fit | Residual |
| (g/210L) | (g/210L) | (g/210L) | (g/210L) | (g/210L) | (g/210L) |
| 0.000 | 0.000 | -0.0003 | 0.000 | 0.000 | -0.0004 |
| 0.040 | 0.040 | 0.0004 | 0.040 | 0.039 | 0.0006 |
| 0.080 | 0.079 | 0.0008 | 0.080 | 0.079 | 0.0006 |
| 0.100 | 0.101 | -0.0009 | 0.100 | 0.101 | -0.0008 |
| 0.300 | 0.300 | 0.0000 | 0.300 | 0.300 | 0.0001 |
| ----- | | | ----- | | |

| | <<<<< 3um >>>>> | <<<<< 9um >>>>> |
|--|-----------------|-----------------|
| | ----- | ----- |
| Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1 | | |
| Sample | | |
| Sample #1 | 3477.00 | 3371.00 |
| Sample #2 | 3368.00 | 3369.00 |
| Sample #3 | 3353.00 | 3361.00 |
| Sample #4 | 3460.00 | 3384.00 |
| Avg | 3393.6667 | 3371.3333 |
| STD DEV | 57.9339 | 11.6762 |
| REL STD DEV | 1.707 | 0.346 |
| H2O adjust (mg/l*10k) | 416 | 438 |

Atmospheric Pressure = 960

*****CALIBRATION SUCCESSFUL*****

Alu Muel
 05 June 2025