



**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION**

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 5937

A. General Setup and Checks:

1. Initial start up
 - ☒ Power On
 - ☒ 5-minute countdown during initial warm-up
 - ☒ 1-minute countdown to begin test after standby mode
2. ☒ Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 - ☒ Level 1 (Operator)
 - ☒ Level 2 (Field Inspector)
 - ☒ Level 3 (Lab Personnel)
4. ☒ Printer Setup (Level 3,S,P)
5. ☒ Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

1. ☒ Internal Printer Test (Level 1,P). Sign and attach record.
 - ☒ Paper Feed
2. ☒ External Printer Test (Level 1,P). Sign and attach record.
3. Barcode Reader (Level 1,B) (Tone sounds after each)
 - ☒ Operator ID
 - ☒ Invalid ID
 - ☒ ND License
 - ☒ R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

1. ☒ Continuous Air Blank (Level 2,A)
2. ☒ Diagnostic (Level 2,D). Sign and attach test record.
3. ☒ Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
 - Inst Ser#: 80-005937 DSP Ver: 1362.14
 - Core Ver: 8164.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 12965 s: 0 C: 47.00
9: 13402 s: 2 B: 45

☒ DVM (3 & 9) between 9,000-16,000 counts

☒ RFI indicator with keyed radio: 0 → 1

2. Function (F) Keys:

☒ F5 – Tone

☒ F9 - Fan

☒ F6 – Green LED

☒ F10 – Pump

☒ F7 – Red LED

☒ F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 955 CF: 1.05

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 215 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 44.0 I: 0.0

B: 45.0 V: 0.0

☒ Cell Temperature (C) between 47 ± 0.02 °C

☒ Breath Hose Temperature (B) between 35.5 – 47 °C

☒ Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 464 psi Regulator: 500 psi

Display and Regulator ± 50 psi of each other ☒ Yes or No

Gas tank tare necessary? Yes or ☒ No

If Yes, display readings after tare (Level 3,M,C,G):

Display: _____ psi Regulator: _____ psi

7. ☒ General Setup (Level 3,S,G)

8. ☒ Configure Start Test (Level 3,S,T)

9. ☒ Network (Level 3,S,N)

E. Tests (Sign and attach all test records):

1. ☒ Configure wet bath simulator for the following test (Level 1,S)

Zero AC-Wet Calibration Check (Level 1,C):

Sim ser #: MP3003 Lot #: N/A Exp date: N/A

Known value: 0.000 AC

☒ Results are Zero AC

2. ☒ Configure wet bath simulator for the following test (Level 1,S)

Low AC-Wet Calibration Check (Level 1,C):

Known Value ≤ 0.03 AC: 0.020 AC

Sim ser #: MP5321 Lot #: 202311D Exp date: 11/28/25


☒ Results ± 0.005 of known AC

3. ☒ Configure wet bath simulator for the following test (Level 1,S)
 High AC-Wet Calibration Check (Level 1,C):
 Known Value ≥ 0.25 AC: 0.300 AC *aen 6/5/25*
 Sim ser #: MP6035 Lot #: 2023110 Exp date: 11/28/25 *aen 6/5/25*
☒ Results $\pm 5\%$ of known AC 202402C 2114126
4. ☒ Configure dry gas standard for the following tests (Level 1,S)
 Known Value: 0.080 AC
 Gas Cylinder Lot #: 14323080A1 Cylinder #: 68 Exp. Date: 6/5/25
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):
 Sim. Ser #: DRT352 Lot #: ICS8 Exp. Date: None
☒ Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B):
 Sim. Ser #: DR7111 Lot #: ACE4 Exp. Date: None
☒ Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):
 Sim. Ser #: DR7345 Lot #: MET4 Exp. Date: None
☒ Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):
 Sim. Ser #: DRT350 Lot #: IPA4 Exp. Date: None
☒ Display Reads "Interferent Detect"
9. ☒ Ambient Fail (Level 1,B): Display read and test record correct.
10. ☒ Invalid Sample (Level 1,B): Display read and test record correct.
11. ☒ Sequence Aborted (Level 1,B): Display read and test record correct.
12. ☒ ABA (Level 1,B): Test record correct.
13. ☒ Calibration check out of tolerance (CMS mode): Display read and test record correct.
14. ☒ RFI Detected (B,C, or CMS mode): Display read and test record correct.
15. ☒ Tank Pressure Below Minimum (CMS mode): Display reads correct.
16. ☒ Deficient Samples (CMS mode): Display read and test record correct.
17. ☒ Improper Sample (CMS mode): Display read and test record correct.
18. ☒ Difference Too Great (CMS mode): Display read and test record correct.
19. ☒ CMS (CMS mode): Test record correct.
20. Stability Test (Level 3,T,S):
 Ave: 0.0784 Std Dev: 0.0005 Rel Std Dev: 0.65%

Maintenance/Repairs/Notes: NIA

Instrument is acceptable to be used in the field.
If No, state reason(s) why:

☒ Yes or No


Inspector Signature

05 June 2025
Date


Reviewer

05 June 2025
Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 11:00

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMN

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNOPQR

24 characters per line
abcdefghijklmnopqrstuvwx
ABCDEFGHIJKLMNOPQRSTU

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNOPQRSTUWXY123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNQRSTUWXY1234567890ABCDEF


Current Instrument Setup
Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 01524080A4
Standard Cyl #? 21
Standard Expiration? 03/05/2026
Oper No? 133237

Flow Cal. Date: 11/17/2017
Slope 606
Intercept -534563

IR Calibration Date: 11/17/2017
3um 9um

0th Coef(*100): -30534 -23772
1st Coef(*100): 294336 142357
2nd Coef(*100): -906 850
H2O adj(mg/l*10K): 205 250

***** Test End *****


Operator Signature
ANNA NAREHOOD

Remarks:

Internal Print
Test
Form 106-18000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 12:07

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%^&*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 01524080A4
Standard Cyl #? 21
Standard Expiration? 03/05/2026
Oper No? 133237

Flow Cal. Date: 11/17/2017
Slope 606
Intercept -534563

IR Calibration Date: 11/17/2017
 3um 9um

0th Coef(*100):	-30534	-23772
1st Coef(*100):	294336	142357
2nd Coef(*100):	-906	850
H2O adj(mg/l*10k):	205	250

***** Printer Test End *****



Operator Signature
ANNA NAREHOOD

Remarks: External print Test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 12:10

Test	AC	Time
01 Room Air	0.000	12:11
02 *Subject Test	REF*	12:11
03 Room Air	0.000	12:12

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = REFUSAL
Dr. Lic. = ND/DIS821456
Lot No = 01524080A4
Cyl No = 21
Expiration Date = 03/05/2026
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Refusal

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 12:13

DIAGNOSTIC

Voltage/Current Test	Pass
RAM Test	Pass
EEPROM Checksum Test	Pass
Real Time Clock Test	Pass
DSP Test	Pass
Analytical Stability Test	Pass
Modem Test	Pass
Temp Regulation Test	Pass

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Diagnostic

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 12:39

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:40
02 Std. Sol.	0.000	12:41
03 Room Air	0.000	12:41
04 Std. Sol.	0.000	12:42
05 Room Air	0.000	12:43
06 Std. Sol.	0.000	12:43
07 Room Air	0.000	12:44


08 Sim Temp = 34.0°C

Simul Ser No = MP3003

Std Sol No = NA

County = 08

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Zero AC Check -0.000AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 12:48

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:49
02 Std. Sol.	0.020	12:52
03 Room Air	0.000	12:52
04 Std. Sol.	0.019	12:53
05 Room Air	0.000	12:54
06 Std. Sol.	0.020	12:54
07 Room Air	0.000	12:55

08 Sim Temp = 34.0°C

Simul Ser No = MP5321
Std Sol No = 202311D
County = 08 Oper No. = 133237


Operator Signature
ANNA NAREHOOD

Remarks: LOW AC Check - 0.020AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 12:57

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:58
02 Std. Sol.	0.298	12:59
03 Room Air	0.000	12:59
04 Std. Sol.	0.301	13:00
05 Room Air	0.000	13:01
06 Std. Sol.	0.301	13:01
07 Room Air	0.000	13:02

08 Sim Temp = 34.0°C

Simul Ser No = MP6035

Std Sol No = ~~202311D~~ 202402C

County = 08

Oper No. = 133237

aen 05 June 2025



Operator Signature
ANNA NAREHOOD

Remarks: *High AC check - 0.300AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:06

Test	AC	Time
01 Room Air	0.000	13:08
02 *Subject Test	INT*	13:08
03 Room Air	0.000	13:09

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT ICS
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10AL + 0.05% Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:09

Test	AC	Time
01 Room Air	0.000	13:10
02 *Subject Test	INT*	13:11
03 Room Air	0.000	13:12

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT ACE
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.05% Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:12

Test	AC	Time
01 Room Air	0.000	13:13
02 *Subject Test	INT*	13:14
03 Room Air	0.000	13:14

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT MET
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10% Methanol

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:15

Test	AC	Time
01 Room Air	0.000	13:17
02 *Subject Test	INT*	13:17
03 Room Air	0.000	13:18

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT IPA
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237


Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10% Isoopropanol

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:18

Test	AC	Time
01 Room Air	AMB*	13:19
02 Room Air	0.000	13:20

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = AMBIENT
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Ambient Fail

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:21

Test	AC	Time
01 Room Air	0.000	13:22
02 Invalid Sample	X.XXX	13:22
03 Room Air	0.000	13:23

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INVALID
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Invalid Sample

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:23

Test	AC	Time
01 Room Air	ABT*	13:24
02 Room Air	0.000	13:24

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = ABORT SEQUENCE
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Sequence Aborted.

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:30

Test	AC	Time
01 Room Air	0.000	13:31
02 Subject Test 1	0.000	13:32
03 Room Air	0.000	13:32
04 Reported AC	0.000	13:32

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = ABA
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: ABA TEST

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:33

Test	AC	Time
01 Diagnostic	OK	13:33
02 Room Air	0.000	13:34
03 Subject Test 1	0.000	13:34
04 Room Air	0.000	13:38
05 Std. Gas	0.000*	13:39
06 Room Air	0.000	13:39

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = CAL CHECK OOT
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Calibration Check out
of tolerance

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:40
02 Std. Gas	RFI*	13:40
03 Room Air	0.000	13:41

*Invalid Test
Inhibited - RFI

Lot No = 14323080A1
Cyl No = 68
Exp Date = 06/05/2025
County = 08 Oper No. = 133237


Operator Signature
ANNA NAREHOOD

Remarks: RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 13:41

Test	AC	Time
01 Diagnostic	OK	13:42
02 Room Air	0.000	13:43
03 *Subject Test	0.000*	13:46
04 Room Air	0.000	13:48
05 Std. Gas	0.078	13:50
06 Room Air	0.000	13:51
07 *Subject Test	0.000*	13:54
08 Room Air	0.000	13:55

No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = DEFICIENT
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Deficient Samples

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 14:01

Test	AC	Time
01 Diagnostic	OK	14:02
02 Room Air	IPS*	14:02
03 Room Air	0.000	14:03

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = IMPROPER
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Improper Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 14:04

Test	AC	Time
01 Diagnostic	OK	14:05
02 Room Air	0.000	14:05
03 Subject Test 1	0.000*	14:06
04 Room Air	0.000	14:08
05 Std. Gas	0.079	14:09
06 Room Air	0.000	14:11
07 Subject Test 2	0.103*	14:11
08 Room Air	0.000	14:12

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.


Operator Signature
ANNA NAREHOOD

Remarks: Difference Too Great

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 14:13

Test	AC	Time
01 Diagnostic	OK	14:14
02 Room Air	0.000	14:14
03 Subject Test 1	0.000	14:15
04 Room Air	0.000	14:18
05 Std. Gas	0.078	14:19
06 Room Air	0.000	14:21
07 Subject Test 2	0.000	14:21
08 Room Air	0.000	14:22
09 Reported AC	0.000	14:15

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = CMS TEST
Dr. Lic. = ND/DIS821456
Lot No = 14323080A1
Cyl No = 68
Expiration Date = 06/05/2025
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD


Remarks: CMS Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = TOXL 8164.14.00 09/16
06/05/2025 14:23

Test	AC	Time
01 Room Air	0.000	14:24
02 Std. Gas	0.078	14:24
03 Room Air	0.000	14:25
04 Std. Gas	0.078	14:25
05 Room Air	0.000	14:26
06 Std. Gas	0.078	14:26
07 Room Air	0.000	14:26
08 Std. Gas	0.078	14:27
09 Room Air	0.000	14:27
10 Std. Gas	0.079	14:28
11 Room Air	0.000	14:28
12 Std. Gas	0.079	14:29
13 Room Air	0.000	14:29
14 Std. Gas	0.078	14:30
15 Room Air	0.000	14:30
16 Std. Gas	0.079	14:31
17 Room Air	0.000	14:31
18 Std. Gas	0.078	14:31
19 Room Air	0.000	14:32
20 Std. Gas	0.079	14:32
21 Room Air	0.000	14:33
22 Std. Gas	0.078	14:33
23 Room Air	0.000	14:34
24 Std. Gas	0.078	14:34
25 Room Air	0.000	14:35
26 Std. Gas	0.079	14:35
27 Room Air	0.000	14:36
28 Std. Gas	0.078	14:36
29 Room Air	0.000	14:37
30 Std. Gas	0.079	14:37
31 Room Air	0.000	14:38
32 Std. Gas	0.078	14:38
33 Room Air	0.000	14:39
34 Std. Gas	0.078	14:39
35 Room Air	0.000	14:39
36 Std. Gas	0.079	14:40
37 Room Air	0.000	14:40
38 Std. Gas	0.078	14:41
39 Room Air	0.000	14:41
40 Std. Gas	0.079	14:42
41 Room Air	0.000	14:42
42 Cal Check Stats:		
43 Average	0.0784	
44 Std Dev	0.0005	
45 Rel Std Dev	0.65%	


05 June 2025

Operator Signature
ANNA NAREHOOD

Remarks: Stability Test

Anna Narehood
05 June 2025

Form 106-I8000