	NORTH DAKOTA OFFICE OF ATTORNEY GE CRIME LABORATORY DIVISION	NERAL
	INTOXILYZER® 8000 ANNUAL INSPECT	ΓΙΟΝ
Intoxilyze	r [®] 8000 Serial Number: 80-00_ <u>5951</u> Inspe	ction Location: TOXL
 A. Pre 1. 2. 3. 4. 5. 6. 	e-Inspection Items with Instrument: Gas Cylinder Yes of No (If Yes, Lot # Keys Yes of No Power Cord Yes or No Download Data Upload Operator File Current Location Code: <u>TOXL</u> Battery Check Was the external battery pack replaced? O-Rings Replaced Simulator O-Ring Replaced Breath Tube O-Ring	Cyl. #) Yes or No Yes or No Yes or No
B. Ge 1. 2. 3. 4. 5.	 Anternal Setup and Checks: Diagnostics passed and instrument in "Ready" modeling Breath tube heated Date, time and location code (Level 2,E). Re-set in Time Zone: CST or MDT (Time on test record) Print test (Level 1,P). Sign and attach test record. Print test (Level 1,P). Sign and attach test record. Tank monitor (Level 3,D,G). Display: <u>482</u> psi Regulator: <u>500</u> Display and Regulator ± 50 psi of each other Gas tank tare necessary? If Yes, display readings after tare (Level 3,M,C Display: psi Regulator: 	f necessary. Is will be in time zone circled). psi es or No Yes or No 2,G):
1. Intoxil	ests (Sign and attach test records): X Configure simulator for the following test (Level 1 Wet Calibration Check - Low AC (Level 1,C) Known Value ≤ 0.03 AC: <u>0.020</u> AC Sim. Ser #: <u>MP5321</u> Lot #: <u>2023110</u> Exp. Date: <u>1128</u> X Results ± 0.005 of known AC yzer 8000 Annual Inspection atory Unit: Toxicology Unit - Breath Alcohol Section	,S). Document ID: 11698 Revision: 3 Status: Published
Appro	ved By: Laboratory Director	Date Approved: 03/20/2025

UNCONTROLLED WHEN PRINTED

Page 1 of 2

	2.	Configure simulator for the following test (Wet Calibration Check - High AC (Leve Known Value ≥ 0.25 AC: 0.300 Sim. Ser #: MP6035 Lot #: 202402C Exp. Date: 2 14 26	el 1,C) _AC
	3.	Results ± 5% AC of known AC Configure dry gas standard for the remaini	na tests (Level 1.S).
		Known Value : 0,080 AC	
		Gas Cylinder Lot #:	3
		Cylinder #: 41 Exp. Date: 11 5125	
	4.	X Interferent Check (Level 1,B)	
	4.	Known Value: 0.10 AC + 0.05% Ac	etone
		Sim. Ser #:	
		Lot #: 1058	
		Exp. Date: None	
		Display reads "Interferent Detect"	
	5.	RFI Check (CMS Mode)	
		Display reads "RFI Detect"	
	6.	Dry Calibration Check (Level 1,C)	
		Test 1 0.079 Test 4 0.079 Test 2 0.079 Test 5 0.019	Test 7 0.079
		Test 3 0.079 Test 6 0.019	Test 8 0.079 Test 9 0.079
		Average 0.010	
		Results ± 0.005 AC of known AC	
		X. cooline 2 choose i to children i to	
D.	Rem	narks/Maintenance: NIA	
<u></u>			
Insti	rument	is acceptable to be used in the field. Yes or No	
If No	o, state	reason(s) why:	
		())	
<u></u>		2	
lf Y€	es, chai	nge location code back to A.4 unless A.4 is TOX	
		IN (MILLO	18 June 2025
Insp	ector S	Signature	Date
On		c Pritschuler	18 June 2025 Date 20 June 2025
		C THIS CHARTER	Date
rev	iewer		
		er 8000 Annual Inspection ory Unit: Toxicology Unit - Breath Alcohol Section	Document ID: 11698 Revision: 3 Status: Published
	Laborato	by one reaction by one broad motion of soliding	orardo, i abilorida

Date Approved: 03/20/2025

Approved By: Laboratory Director

UNCONTROLLED WHEN PRINTED

Page 2 of 2

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005951 Location = TOXL 8164.14.00 09/16 06/18/2025 10:01 abcdefqhijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopgrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Wet Standard Value? 0.400 Standard Solution #? 202312D Sim Serial #? MP5318 Oper No? 133237 Flow Cal. Date: 06/21/2013 Slope 693 Intercept -721277 IR Calibration Date: 04/03/2024 3um 9um _____ 0th Coef(*100):-207491st Coef(*100):2625992nd Coef(*100):2775 -22141135404 1304 H2O adj(mg/l*10k): 463445 ********** Printer Test End **********

Operator Signature ANNA NAREHOOD

Remarks:

Test

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005951Location = TOXL8164.14.00 09/1606/18/202510:17

		Ĺ	WET CAL CHECK		
Test			AC	Time	
01	Room	Air	0.000	10:18	
02	Std.	Sol.	0.020	10:19	
03	Room	Air	0.000	10:19	
04	Std.	Sol.	0.020	10:20	
05	Room	Air	0.000	10:21	
06	Std.	Sol.	0.020	10:21	
07	Room	Air	0.000	10:22	

08 Sim Temp = 34.0°C

Simul Ser No = 1	MP5321
Std Sol No = 202	2311D
County = 08	Oper No. = 133237

Operator Signature ANNA NAREHOOD

LOW AC Check - 0.020AC Remarks:

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005951
Location = TOXL	8164.14.00 09/16
06/18/2025	10:23

			WET CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	10:24
02	Std.	Sol.	0.290	10:25
03	Room	Air	0.000	10:26
04	Std.	Sol.	0.292	10:26
05	Room	Air	0.000	10:27
06	Std.	Sol.	0.292	10:28
07	Room	Air	0.000	10:28

08 Sim Temp = 34.0°C

Simul Ser No	= MP6035			
Std Sol No =	202402C			
County = 08		Oper No.	=	133237

Operator Signature ANNA NAREHOOD

Remarks: High AC Check -0.300AC

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005951 Location = TOXL 8164.14.00 09/16 06/18/2025 10:29 AC Time Test 10:30 01 Room Air 0.000 INT* 10:31 02 *Subject Test 03 Room Air 0.000 10:31 *Invalid Test Interferent Detected Sub Name = DISCOVER, THE SPIRIT Sub DOB = 01/01/1982Sub Sex = Female Weight = 150Cit = INTERFERENT ICS Test = OTHDr. Lic. = ND/DIS821456Lot No = 28423080A3Cyl No = 41Expiration Date = 11/05/2025County = 08Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks: Interferent Check 0,10AC + 0.05°10 Acetone

NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-005951 8164.14.00 09/16 Location = TOXL15:32 06/18/2025 Time AC Test 15:33 OK 01 Diagnostic 15:34 02 Room Air 0.000 15:34 RFI* 03 *Subject Test 15:34 04 Room Air 0.000 *Invalid Test Inhibited - RFI Sub Name = DISCOVER, THE SPIRIT Sub DOB = 01/01/1982Sub Sex = Female Weight = 150Cit = RFI CHECK Test = OTHDr. Lic. = ND/DIS821456Lot No = 28423080A3Cyl No = 41Expiration Date = 11/05/2025Oper No. = 133237 County = 08

Intoxilyzer Test Record and Checklist

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature ANNA NAREHOOD

Remarks:

Check

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005951Location = TOXL8164.14.00 09/1606/18/202512:04

			DRY CAL CHECK		
Test			AC	Time	
01	Room	Air	0.000	12:05	
02	Std.	Gas	0.079	12:05	
03	Room	Air	0.000	12:06	
04	Std.	Gas	0.079	12:06	
05	Room	Air	0.000	12:07	
	Std.		0.079	12:07	
07	Room	Air	0.000	12:08	

Lot No = 28423080A3 Cyl No = 41 Exp Date = 11/05/2025 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

y cal check #1-3

Remarks:

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005951Location = TOXL8164.14.00 09/1606/18/202512:10

			DRY CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	12:10
02	Std.	Gas	0.079	12:11
03	Room	Air	0.000	12:11
04	std.	Gas	0.079	12:12
05	Room	Air	0.000	12:12
	Std.		0.079	12:13
	Room		0.000	12:13
0.				

Lot No = 28423080A3 Cyl No = 41 Exp Date = 11/05/2025 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

cal check # 4-6

Remarks: DW

Form 106-18000

1

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005951Location = TOXL8164.14.00 09/1606/18/202512:14

		[DRY CAL CHECK	
Te	est		AC	Time
01	Room	Air	0.000	12:14
02	Std.	Gas	0.079	12:15
03	Room	Air	0.000	12:15
04	Std.	Gas	0.079	12:15
05	Room	Air	0.000	12:16
06	std.	Gas	0.079	12:16
	Room		0.000	12:17

Lot No = 28423080A3 Cyl No = 41 Exp Date = 11/05/2025 County = 08 Oper No. = 133237

Operator Signature ANNA NAREHOOD

Remarks:

calcherk # 7-9