



**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION**

INTOXILYZER® 8000 INITIAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 4943

A. General Setup and Checks:

1. Initial start up
 - ☒ Power On
 - ☒ 5-minute countdown during initial warm-up
 - ☒ 1-minute countdown to begin test after standby mode
2. ☒ Display Intensity (Ctrl + ↑ or ↓)
3. Password Security:
 - ☒ Level 1 (Operator)
 - ☒ Level 2 (Field Inspector)
 - ☒ Level 3 (Lab Personnel)
4. ☒ Printer Setup (Level 3,S,P)
5. ☒ Forms Load Completed (CMI, Inc., CoBrA Software)

B. Level 1 Setup and Checks:

1. ☒ Internal Printer Test (Level 1,P). Sign and attach record.
 - ☒ Paper Feed
2. ☒ External Printer Test (Level 1,P). Sign and attach record.
3. Barcode Reader (Level 1,B) (Tone sounds after each)
 - ☒ Operator ID
 - ☒ Invalid ID
 - ☒ ND License
 - ☒ R key: Refusal. Sign and attach record.

C. Level 2 Setup and Checks:

1. ☒ Continuous Air Blank (Level 2,A)
2. ☒ Diagnostic (Level 2,D). Sign and attach test record.
3. ☒ Prelim Data Entry (Level 2,E)
4. Ser Number/Version Info (Level 2,V):
Inst Ser#: 80-004943 DSP Ver: 1362.14
Core Ver: 8164.14 CPLD Ver: 7

D. Level 3 Setup and Checks:

1. DVM Monitor (Level 3,D,D):

3: 12779 s: 0 C: 47.06
9: 12780 s: 3 B: 44

☒ DVM (3 & 9) between 9,000-16,000 counts

☒ RFI indicator with keyed radio: 0 → 1

2. Function (F) Keys:

☒ F5 – Tone

☒ F9 - Fan

☒ F6 – Green LED

☒ F10 – Pump

☒ F7 – Red LED

☒ F11 – Solenoid

3. Atmospheric Monitor (Level 3,D,A): Atmo Sensor, P: 945 CF: 1.07

4. Flow Monitor (Level 3,D,F): Flow Sensor, R: 222 F: 0.000 L/s

5. Temperature Monitor (Level 3,D,T):

C: 47.0 A: 49.0 I: 0.0

B: 44.9 V: 0.0

☒ Cell Temperature (C) between 47 ± 0.02 °C

☒ Breath Hose Temperature (B) between 35.5 – 47 °C

☒ Internal Ambient Temperature (A) above 41 °C

6. Tank Monitor (Level 3,D,G):

Display: 972 psi Regulator: 1000 psi

Display and Regulator ± 50 psi of each other ☒ Yes or No

Gas tank tare necessary? ☒ Yes or ☒ No

If Yes, display readings after tare (Level 3,M,C,G):

Display: _____ psi Regulator: _____ psi

7. ☒ General Setup (Level 3,S,G)

8. ☒ Configure Start Test (Level 3,S,T)

9. ☒ Network (Level 3,S,N)

E. Tests (Sign and attach all test records):

1. ☒ Configure wet bath simulator for the following test (Level 1,S)

Zero AC-Wet Calibration Check (Level 1,C):

Sim ser #: MP3003 Lot #: NA Exp date: N/A

Known value: 0.000 AC

☒ Results are Zero AC

2. ☒ Configure wet bath simulator for the following test (Level 1,S)

Low AC-Wet Calibration Check (Level 1,C):

Known Value ≤ 0.03 AC: 0.020 AC

Sim ser #: MP5321 Lot #: 20231D Exp date: 11/28/25

☒ Results ± 0.005 of known AC

3. ☒ Configure wet bath simulator for the following test (Level 1,S)
High AC-Wet Calibration Check (Level 1,C):
Known Value ≥ 0.25 AC: 0.300 AC
Sim ser #: MP6035 Lot #: 202402C Exp date: 2/14/26
☒ Results $\pm 5\%$ of known AC
4. ☒ Configure dry gas standard for the following tests (Level 1,S)
Known Value: 0.080 AC
Gas Cylinder Lot #: 28423080A3 Cylinder #: 41 Exp. Date: 11/5/25
5. 0.10% Ethanol + 0.05% Acetone (Level 1,B):
Sim. Ser #: DR7352 Lot #: ICS8 Exp. Date: None
☒ Display reads "Interferent Detect"
6. 0.05% Acetone (Level 1,B): AEN 6/19/25
Sim. Ser #: DR7111 Lot #: AGE4 ACE4 Exp. Date: None
☒ Display reads "Interferent Detect"
7. 0.10% Methanol (Level 1,B):
Sim. Ser #: DR7345 Lot #: MET4 Exp. Date: None
☒ Display Reads "Interferent Detect"
8. 0.10% Isopropanol (Level 1,B):
Sim. Ser #: DR7350 Lot #: IPA4 Exp. Date: None
☒ Display Reads "Interferent Detect"
9. ☒ Ambient Fail (Level 1,B): Display read and test record correct.
10. ☒ Invalid Sample (Level 1,B): Display read and test record correct.
11. ☒ Sequence Aborted (Level 1,B): Display read and test record correct.
12. ☒ ABA (Level 1,B): Test record correct.
13. ☒ Calibration check out of tolerance (CMS mode): Display read and test record correct.
14. ☒ RFI Detected (B,C, or CMS mode): Display read and test record correct.
15. ☒ Tank Pressure Below Minimum (CMS mode): Display reads correct.
16. ☒ Deficient Samples (CMS mode): Display read and test record correct.
17. ☒ Improper Sample (CMS mode): Display read and test record correct.
18. ☒ Difference Too Great (CMS mode): Display read and test record correct.
19. ☒ CMS (CMS mode): Test record correct.
20. Stability Test (Level 3,T,S):
Ave: 0.0800 Std Dev: 0.0005 Rel Std Dev: 0.64%

Maintenance/Repairs/Notes: N/A

Instrument is acceptable to be used in the field.
If No, state reason(s) why:

Yes or No

[Signature]

Inspector Signature

19 June 2025

Date

[Signature]

Reviewer

20 June 2025

Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMi, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 11:25

***** Printer Test *****

12 char/line
abcdefghijkl
ABCDEFGHIJKL

16 char per line
abcdefghijklmnop
ABCDEFGHIJKLMNPO

21 chars per line
abcdefghijklmnopqrstu
ABCDEFGHIJKLMNOPQRSTU

24 characters per line
abcdefghijklmnopqrstuvw
ABCDEFGHIJKLMNOPQRSTUWVX

32 characters per line
abcdefghijklmnopqrstuvwxy123456
ABCDEFGHIJKLMNOPQRSTUWVXY123456

42 characters per line
abcdefghijklmnopqrstuvwxy1234567890abcdef
ABCDEFGHIJKLMNOPQRSTUWVXY1234567890ABCDEF

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On MemFull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 08425080A1
Standard Cyl #? 24
Standard Expiration? 05/05/2027
Oper No? 133237

Flow Cal. Date: 04/06/2020
Slope 679
Intercept -694345

IR Calibration Date: 04/06/2020
3um 9um

| | | |
|--------------------|--------|--------|
| 0th Coef(*100): | -23601 | -21377 |
| 1st Coef(*100): | 255715 | 132279 |
| 2nd Coef(*100): | 1319 | 1274 |
| H2O adj(ng/l*10K): | 424 | 419 |

***** Test End *****



Operator Signature
ANNA NAREHOOD

Remarks:

Internal Print
Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 12:29

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 08425080A1
Standard Cyl #? 24
Standard Expiration? 05/05/2027
Oper No? 133237

Flow Cal. Date: 04/06/2020
Slope 679
Intercept -694345

IR Calibration Date: 04/06/2020
 3um 9um

| | | |
|--------------------|--------|--------|
| 0th Coef(*100): | -23601 | -21377 |
| 1st Coef(*100): | 255715 | 132279 |
| 2nd Coef(*100): | 1319 | 1274 |
| H2O adj(mg/l*10k): | 424 | 419 |

***** Printer Test End *****



Operator Signature
ANNA NAREHOOD

Remarks: External Print
Test

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 12:29

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 12:30 |
| 02 *Subject Test | REF* | 12:30 |
| 03 Room Air | 0.000 | 12:31 |

*Subject Test Refused

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = REFUSAL
Dr. Lic. = ND/DIS821456
Lot No = 08425080A1
Cyl No = 24
Expiration Date = 05/05/2027
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Refusal

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 12:31

DIAGNOSTIC

| | |
|---------------------------|------|
| Voltage/Current Test | Pass |
| RAM Test | Pass |
| EEPROM Checksum Test | Pass |
| Real Time Clock Test | Pass |
| DSP Test | Pass |
| Analytical Stability Test | Pass |
| Modem Test | Pass |
| Temp Regulation Test | Pass |

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *Diagnostic*

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 12:36

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 12:37 |
| 02 Std. Sol. | 0.000 | 12:38 |
| 03 Room Air | 0.000 | 12:38 |
| 04 Std. Sol. | 0.000 | 12:39 |
| 05 Room Air | 0.000 | 12:40 |
| 06 Std. Sol. | 0.000 | 12:40 |
| 07 Room Air | 0.000 | 12:41 |


08 Sim Temp = 34.0°C

Simul Ser No = MP3003

Std Sol No = NA

County = 08

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks:

Zero AC check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 12:41

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 12:43 |
| 02 Std. Sol. | 0.019 | 12:44 |
| 03 Room Air | 0.000 | 12:44 |
| 04 Std. Sol. | 0.019 | 12:45 |
| 05 Room Air | 0.000 | 12:46 |
| 06 Std. Sol. | 0.019 | 12:46 |
| 07 Room Air | 0.000 | 12:47 |


08 Sim Temp = 34.0°C

Simul Ser No = MP5321

Std Sol No = 202311D

County = 08

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks:

Low AC check -0.020AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 12:47

WET CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 12:54 |
| 02 Std. Sol. | 0.296 | 12:55 |
| 03 Room Air | 0.000 | 12:56 |
| 04 Std. Sol. | 0.296 | 12:57 |
| 05 Room Air | 0.000 | 12:57 |
| 06 Std. Sol. | 0.295 | 12:58 |
| 07 Room Air | 0.000 | 12:58 |


08 Sim Temp = 34.0°C

Simul Ser No = MP6035

Std Sol No = 202402C

County = 08

Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: High Ac Check - 0.300AC

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:00

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:01 |
| 02 *Subject Test | INT* | 13:02 |
| 03 Room Air | 0.000 | 13:02 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT ICS
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10AC + 0.050% Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:02

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:03 |
| 02 *Subject Test | INT* | 13:04 |
| 03 Room Air | 0.000 | 13:04 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT ACE
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent Check
0.05% Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:04

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:05 |
| 02 *Subject Test | INT* | 13:05 |
| 03 Room Air | 0.000 | 13:06 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT MET
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Interferent check
0.10% Methanol

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:06

| Test | AC | Time |
|------------------|-------|-------|
| 01 Room Air | 0.000 | 13:07 |
| 02 *Subject Test | INT* | 13:07 |
| 03 Room Air | 0.000 | 13:08 |

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INTERFERENT IPA
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: *interferent check*
0.100% Isopropanol

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:08

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | AMB* | 13:09 |
| 02 Room Air | 0.000 | 13:09 |

*Invalid Test
Check Ambient Conditions

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = AMBIENT
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: Ambient Fail

Form 106-I8000


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:11

| Test | AC | Time |
|-------------------|-------|-------|
| 01 Room Air | 0.000 | 13:11 |
| 02 Invalid Sample | X.XXX | 13:12 |
| 03 Room Air | 0.000 | 13:12 |

*Invalid Test - Mouth Alcohol

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = INVALID
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: invalid sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:13

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | ABT* | 13:13 |
| 02 Room Air | 0.000 | 13:14 |

*Sequence Aborted

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = SEQUENCE ABORT
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks:

Sequence Aborted

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:19

| Test | AC | Time |
|-------------------|-------|-------|
| 01 Room Air | 0.000 | 13:20 |
| 02 Subject Test 1 | 0.000 | 13:20 |
| 03 Room Air | 0.000 | 13:21 |
| 04 Reported AC | 0.000 | 13:20 |

No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = ABA TEST
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237



Operator Signature
ANNA NAREHOOD

Remarks: ABA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:21

| Test | AC | Time |
|-------------------|--------|-------|
| 01 Diagnostic | OK | 13:22 |
| 02 Room Air | 0.000 | 13:23 |
| 03 Subject Test 1 | 0.000 | 13:23 |
| 04 Room Air | 0.000 | 13:25 |
| 05 Std. Gas | 0.000* | 13:26 |
| 06 Room Air | 0.000 | 13:27 |

*Cal Check Out of Tolerance

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = CAL CHECK OOT
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Calibration check out
of Tolerance

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:27

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 13:28 |
| 02 Room Air | RFI* | 13:28 |
| 03 Room Air | 0.000 | 13:29 |

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = RFI CHECK
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:29

| Test | AC | Time |
|------------------|--------|-------|
| 01 Diagnostic | OK | 13:30 |
| 02 Room Air | 0.000 | 13:30 |
| 03 *Subject Test | 0.000* | 13:33 |
| 04 Room Air | 0.000 | 13:36 |
| 05 Std. Gas | 0.079 | 13:37 |
| 06 Room Air | 0.000 | 13:38 |
| 07 *Subject Test | 0.000* | 13:41 |
| 08 Room Air | 0.000 | 13:43 |


No RFI Detected

*Deficient Sample - Value Printed was
Highest Obtained

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = DEFICIENT
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Deficient Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:46

| Test | AC | Time |
|---------------|-------|-------|
| 01 Diagnostic | OK | 13:46 |
| 02 Room Air | IPS* | 13:46 |
| 03 Room Air | 0.000 | 13:47 |

*Improper Sample

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = IMPROPER
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks:

Improper Sample

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:47

| Test | AC | Time |
|-------------------|--------|-------|
| 01 Diagnostic | OK | 13:48 |
| 02 Room Air | 0.000 | 13:49 |
| 03 Subject Test 1 | 0.000* | 13:49 |
| 04 Room Air | 0.000 | 13:51 |
| 05 Std. Gas | 0.080 | 13:52 |
| 06 Room Air | 0.000 | 13:53 |
| 07 Subject Test 2 | 0.077* | 13:54 |
| 08 Room Air | 0.000 | 13:55 |

*Difference Too Great

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = DIFF TOO GREAT
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: Difference Too Great

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 13:56

| Test | AC | Time |
|-------------------|-------|-------|
| 01 Diagnostic | OK | 13:57 |
| 02 Room Air | 0.000 | 13:57 |
| 03 Subject Test 1 | 0.000 | 13:57 |
| 04 Room Air | 0.000 | 13:59 |
| 05 Std. Gas | 0.079 | 14:01 |
| 06 Room Air | 0.000 | 14:02 |
| 07 Subject Test 2 | 0.000 | 14:03 |
| 08 Room Air | 0.000 | 14:04 |
| 09 Reported AC | 0.000 | 13:57 |

Difference OK
No RFI Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 01/01/1982
Sub Sex = Female Weight = 150
Test = OTH Cit = CMS TEST
Dr. Lic. = ND/DIS821456
Lot No = 28423080A3
Cyl No = 41
Expiration Date = 11/05/2025
County = 08 Oper No. = 133237

20 minute waiting period ascertained? Y
Clean Mouthpiece used and disposed? Y
Clean Mouthpiece used and disposed? Y

I followed the Approved Method and the
instructions displayed by the Intoxilyzer
in conducting this test.



Operator Signature
ANNA NAREHOOD

Remarks: CMS Test

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004943
Location = TOXL 8164.14.00 09/16
06/19/2025 14:06

| Test | AC | Time |
|---------------------|--------|-------|
| 01 Room Air | 0.000 | 14:06 |
| 02 Std. Gas | 0.079 | 14:07 |
| 03 Room Air | 0.000 | 14:07 |
| 04 Std. Gas | 0.080 | 14:08 |
| 05 Room Air | 0.000 | 14:08 |
| 06 Std. Gas | 0.079 | 14:08 |
| 07 Room Air | 0.000 | 14:09 |
| 08 Std. Gas | 0.080 | 14:09 |
| 09 Room Air | 0.000 | 14:10 |
| 10 Std. Gas | 0.080 | 14:10 |
| 11 Room Air | 0.000 | 14:11 |
| 12 Std. Gas | 0.080 | 14:11 |
| 13 Room Air | 0.000 | 14:11 |
| 14 Std. Gas | 0.080 | 14:12 |
| 15 Room Air | 0.000 | 14:12 |
| 16 Std. Gas | 0.080 | 14:13 |
| 17 Room Air | 0.000 | 14:13 |
| 18 Std. Gas | 0.080 | 14:14 |
| 19 Room Air | 0.000 | 14:14 |
| 20 Std. Gas | 0.080 | 14:14 |
| 21 Room Air | 0.000 | 14:15 |
| 22 Std. Gas | 0.080 | 14:15 |
| 23 Room Air | 0.000 | 14:16 |
| 24 Std. Gas | 0.080 | 14:16 |
| 25 Room Air | 0.000 | 14:17 |
| 26 Std. Gas | 0.080 | 14:17 |
| 27 Room Air | 0.000 | 14:17 |
| 28 Std. Gas | 0.079 | 14:18 |
| 29 Room Air | 0.000 | 14:18 |
| 30 Std. Gas | 0.080 | 14:19 |
| 31 Room Air | 0.000 | 14:19 |
| 32 Std. Gas | 0.080 | 14:20 |
| 33 Room Air | 0.000 | 14:20 |
| 34 Std. Gas | 0.081 | 14:20 |
| 35 Room Air | 0.000 | 14:21 |
| 36 Std. Gas | 0.080 | 14:21 |
| 37 Room Air | 0.000 | 14:22 |
| 38 Std. Gas | 0.080 | 14:22 |
| 39 Room Air | 0.000 | 14:23 |
| 40 Std. Gas | 0.081 | 14:23 |
| 41 Room Air | 0.000 | 14:23 |
| 42 Cal Check Stats: | | |
| 43 Average | 0.0800 | |
| 44 Std Dev | 0.0005 | |
| 45 Rel Std Dev | 0.64% | |

Operator Signature
ANNA NAREHOOD

Remarks:

Stability check

A handwritten signature in blue ink, appearing to be 'Anna Narehood', written in a cursive style.

Form 106-I8000