

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Social Number					
Serial Number 80-004196 Reason for Install/Repair	Instrument Location 1+e++	(Adams (ounty)			
N Install Ages B	r Location Change				
Check When Done:					
1. Surge Protector Installed/Property Grounded.					
Z. Telephone Line Connected to Intoxilyzer® 8000.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
€ 6. Run Tests:					
A. Print Test (Level1, Function P).	A. Print Test (Level1, Function P).				
B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).					
7. Repair and/or Maintenance Performed (if any):	Inspection				
8. Complete the Top Portion of the Intoxilyzer® Record (SFN504		ne Intoxilyzer® for Use.			
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.					
10. Send the Following to the Crime Laboratory: Completed Intoxi 104-G), Print Test, ACA Test, and RFI Test.	ilyzer® 8000 Installation and Repa	ir Checkout (SFN59281, Form			
Field Inspector Signature		Date 6-27-25			
Crime Laboratory Use Only					
This installation has been reviewed and the instrument is approved to be u the date the Field Inspector performed the installation. This record on file a County of Burleigh, North Dakota, is certified to be a true and correct copy	at the Office of Attorney General, C	etermine alcohol concentration from Crime Laboratory Division, in the			
Reviewed/Certified By Anna Narehood	nolo	Certified Date			

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-004196 Location = HETT 8164.14.00 09/16 06/27/2025 11:29

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz12345t7890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 08425080A1 14 Standard Expiration? 05/05/2027 Oper No? 133194

Flow Cal. Date: 04/15/2020 Slope 657 Intercept -683120

IR Calibration Date: 04/16/2020 3um 0th Coef(*100): -35275 -26830 1st Coef(*100): 271852 134658 2nd Coef(*100): 2/105. 849 H2O adj(mg/l*10k): 413

******* Printer Test End ********

302

perator Signature JONATHAN PIEHL

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HETT
06/27/2025

SN 80-004196 8164.14.00 09/16 11:37

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	11:39
02	Std.	Gas	0.077	11:39
03	Room	Air	0.000	11:39
04	Std.	Gas	0.077	11:40
05	Room	Air	0.000	11:40
06	Std.	Gas	0.077	11:41
07	Room	Air	0.000	11:41

Lot No = 08425080A1

Cyl No = 14

Exp Date = 05/05/2027

County = 01

Oper No. = 133194

Operator Signature JONATHAN PIEHL

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = HETT 06/27/2025

SN 80-004196 8164.14.00 09/16 11:44

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	11:44
02 Room Air	0.000	11:45

*Invalid Test Inhibited - RFI

Lot No = 08425080A1

Cyl No = 14

Exp Date = 05/05/2027

County = 01

Oper No. = 133194

perator Signature JONATHAN PIEHL

Remarks:

Form 106-I8000