

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 00650 (	Instrument Location Minot AFB		
Reason for Install/Repair			
Install After Receiving From Crime Laboratory	er Location Change		
Other (Specify)			
Check When Done:			
1. Surge Protector Installed/Property Grounded.			
→ 2. Telephone Line Connected to Intoxilyzer® 8000.			
3. Breath Tube Heated.			
又 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).		
★ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
★ 6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).			
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5)	0496, Form 120-G) and Place it by the Ir	ntoxilyzer® for Use.	
	i) at the Intoxilyzer® Location at the Age	ncy.	
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature	N. C.	Date	
Jefler IT Salwy		1 July 2025	
Crime Laboratory Use Only			
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on fi County of Burleigh, North Dakota, is certified to be a true and correct co	le at the Office of Attorney General, Crin	rmine alcohol concentration from ne Laboratory Division, in the	
Reviewed/Certified By Anna Narehood	1 Male	Certified Date 15July2025	

FORM 104-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-006501
Location = MNAF	8164.14.00 09/16
07/01/2025	15:08

\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*()\_+?

Current	Instrument	Setup
---------	------------	-------

Data Entry Mode:	Enabled
Start Test Sequence:	DABACABA
Display Prelim Rslt?	Yes
Display Third Digit?	Yes
Inhib Printer(Y/N)?	No
Display Volume?	No
Disable On Memfull?	Yes
# of Print Copies?	1
Select Std (D/W/I)?	Dry
Standard Value?	0.080
Standard Lot #?	25024080A3
Standard Cyl #?	50
Standard Expiration?	10/05/2026
Oper No?	134040
and the second s	

Flow Cal. Date:	06/13/2019
Slope	651
Intercept	-777821

IR Calibration Date	: 06/13/2019	
	3um	9um
0th Coef(*100):	-20531	-14906
1st Coef(*100):	270252	139749
2nd Coef(*100):	1917	1284
H2O adj(mg/l*10k):	393	321

\*\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*\*

Operator Signature
TYLER MOLLOY

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006501 Location = MNAF 07/01/2025

8164.14.00 09/16 15:00

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	15:04
02	Std.	Gas	0.080	15:05
03	Room	Air	0.000	15:05
04	Std.	Gas	0.079	15:06
05	Room	Air	0.000	15:06
06	Std.	Gas	0.079	15:07
07	Room	Air	0.000	15:07

Lot No = 25024080A3

Cyl No = 50

Exp Date = 10/05/2026

County = 64

Oper No. = 134040

Óperator Signature TYLER MOLLOY

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MNAF
07/01/2025

Alcohol Analyzer SN 80-006501 8164.14.00 09/16 15:10

## DRY CAL CHECK

Test			AC	Time
01	Room	Air	RFI*	15:11
02	Room	Air	RFI*	15:11

\*Invalid Test Inhibited - RFI

Lot No = 25024080A3

Cyl No = 50

Exp Date = 10/05/2026

County = 64

Oper No. = 134040

Operator Signature
TYLER MOLLOY

Remarks:

Form 106-I8000