

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

| Serial Number Instrument Location |
|--|
| Reason for Install/Repair |
| |
| Install After Receiving From Crime Laboratory Install After Location Change |
| Other (Specify) Lostall after repairs |
| Check When Done: |
| 1. Surge Protector Installed/Property Grounded. |
| 2. Telephone Line Connected to Intoxilyzer® 8000. |
| 3. Breath Tube Heated. |
| 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E). |
| 5. Scan/Enter Gas Cylinder Information (Level 1, Function S). |
| 6. Run Tests: |
| A. Print Test (Level1, Function P). |
| B. ACA Test (Level 1, Function C). |
| C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test). |
| 7. Repair and/or Maintenance Performed (if any): |
| 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use. |
| 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency. |
| 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test. |
| Field Inspector Signature Date OTIVA 2005 |
| Crime Laboratory Use Only |
| This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received. |
| Reviewed/Certified By Anna Narehood Certified Date 17July2025 |

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005955
Location = WILL 8164.14.00 09/16
07/16/2025 11:21

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 25024080A3 Standard Cyl #? 56 Standard Expiration? 10/05/2026 Oper No? 132166

Flow Cal. Date: 04/19/2022 Slope 2446 Intercept -3112812

IR Calibration Date: 06/21/2023
3um 9um

Oth Coef(*100): -12629 -18041
1st Coef(*100): 273918 134302
2nd Coef(*100): 3542 1308
H2O adj(mg/l*10k): 397 389

****** Printer Test End *******

Operator Signatur ALYSON MOSES

Remarks:

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WILL
07/16/2025

Alcohol Analyzer SN 80-005955 8164.14.00 09/16 11:23

DRY CAL CHECK

| Τe | est | | AC | Time |
|----|------|-----|-------|-------|
| 01 | Room | Air | 0.000 | 11:24 |
| 02 | Std. | Gas | 0.081 | 11:24 |
| 03 | Room | Air | 0.000 | 11:25 |
| 04 | Std. | Gas | 0.080 | 11:25 |
| 05 | Room | Air | 0.000 | 11:26 |
| 06 | Std. | Gas | 0.080 | 11:26 |
| 07 | Room | Air | 0.000 | 11:27 |

Lot No = 25024080A3

Cyl No = 56

Exp Date = 10/05/2026

County = 53

Oper No. = 132166

Remarks:

Form 106-I8000

Operator Signature
ALYSON MOSES

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 Location = WILL

SN 80-005955 8164.14.00 09/16 11:22

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | RFI* | 11:23 |
| 02 Room Air | 0.000 | 11:23 |

*Invalid Test Inhibited - RFI

Lot No = 25024080A3

Cyl No = 56

07/16/2025

Exp Date = 10/05/2026

County = 53

Oper No. = 132166

Operator Signature ALYSON MOSES

Remarks:

Form 106-I8000