

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-006500 Instrument Location
Reason for Install/Repair
Install After Receiving From Crime Laboratory Install After Location Change
Other (Specify)
Check When Done:
7. Surge Protector Installed/Property Grounded.
2. Telephone Line Connected to Intoxilyzer® 8000.
3. Breath Tube Heated.
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
6. Run Tests:
A. Print Test (Level1, Function P).
B. ACA Test (Level 1, Function C).
Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any):
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.
Field Inspector Signature Date 8/26/25
Crime Laboratory Use Only
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.
Reviewed/Certified By Anna Narehood Certified Date 27August2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006500
Location = HARV 8164.14.00 09/16
08/26/2025 10:53

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? Standard Cyl #? 01524080A4 25 Standard Expiration? 03/05/2026 Oper No? 132537

Flow Cal. Date: 08/18/2015 Slope 681 Intercept -680985

IR Calibration Date: 08/18/2015
3um 9um

Oth Coef(*100): -3148 -17259
1st Coef(*100): 254475 138520
2nd Coef(*100): 1954 1329
H2O adj(mg/l*10k): 193 414

****** Printer Test End *******

Operator Signature GABRIELA URZOLA _8/26/2S

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HARV
08/26/2025

Alcohol Analyzer SN 80-006500 8164.14.00 09/16 11:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:03
02 Std. Gas	0.081	11:04
03 Room Air	0.000	11:04
04 Std. Gas	0.081	11:05
05 Room Air	0.000	11:05
06 Std. Gas	0.081	11:06
07 Room Air	0.000	11:06

Lot No = 01524080A4

Cyl No = 25

Exp Date = 03/05/2026

County = 52

Oper No. = 132537

Operator Signature
GABRIELA URZOLA

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HARV
08/26/2025

Alcohol Analyzer SN 80-006500 8164.14.00 09/16 10:58

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:59
02 Std. Gas	0.081	10:59
03 Room Air	0.000	11:00
04 Std. Gas	0.080	11:00
05 Room Air	0.000	11:01
06 Std. Gas	0.081	11:01
07 Room Air	RFI*	11:01
08 Room Air	RFI*	11:01

*Invalid Test Inhibited - RFI

Lot No = 01524080A4

Cyl No = 25

Exp Date = 03/05/2026

County = 52

Oper No. = 132537

Cable 2 V/2. 2 8/26/29
Operator Signature
GABRIELA URZLA

Remarks:

Form 106-I8000