

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-006506	Instrument Location
Reason for Install/Repair	Morton County Sheriff's Office Mande
Install After Receiving From Crime Laboratory Install	After Location Change
Other (Specify)	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Lo	ocation; Level 2, Function E).
5. Scan/Enter Gas Cylinder Information (Level 1, Function	S).
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mod	e or Level 1, Function B or C; Key Radio During Test).
7. Repair and/or Maintenance Performed (if any): Nor	<i>1€</i>
8. Complete the Top Portion of the Intoxilyzer® Record (SF	N50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 12	20-G) at the Intoxilyzer® Location at the Agency.
10. Send the Following to the Crime Laboratory: Completed 104-G), Print Test, ACA Test, and RFI Test.	Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form
Field Inspector Signature	Date 8-13-2015
Crime Laboratory Use Only	
	to be used for the analysis of breath to determine alcohol concentration from on file at the Office of Attorney General, Crime Laboratory Division, in the t copy of the documents received.
Reviewed/Certified By	Certified Date

Anna Narehood

Certified Date

03Sept2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006506
Location = MAND 8164.14.00 09/16
08/23/2025 05:07

********* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? NoDisable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 25024080A3 Standard Lot #? Standard Cyl #? Standard Expiration? 10/05/2026 133572 Oper No?

Flow Cal. Date: 08/20/2015 Slope 657 Intercept -627881

IR Calibration Date: 08/20/2015

Oth Coef(*100): -22470 -18919
1st Coef(*100): 266505 136109
2nd Coef(*100): 3702 1644
H2O adj(mq/l*10k): 709 510

3um

9um

****** Printer Test End *******

Operator Signature LOGAN MOLLMAN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MAND
08/23/2025

Alcohol Analyzer SN 80-006506 8164.14.00 09/16 05:29

DRY CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	05:30
02	Std.	Gas	0.080	05:30
03	Room	Air	0.000	05:31
04	Std.	Gas	0.080	05:31
05	Room	Air	0.000	05:31
06	Std.	Gas	0.080	05:32
07	Room	Air	0.000	05:32

Lot No = 25024080A3

Cyl No = 10

Exp Date = 10/05/2026

County = 30

Oper No. = 133572

Operator Signature LOGAN MOLLMAN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = MAND
08/23/2025

Alcohol Analyzer SN 80-006506 8164.14.00 09/16 05:33

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	05:34
02 Std. Gas	0.080	05:34
03 Room Air	RFI*	05:34
04 Room Air	0.000	05:35

*Invalid Test Inhibited - RFI

Lot No = 25024080A3

Cyl No = 10

Exp Date = 10/05/2026

County = 30

Oper No. = 133572

Operator Signature LOGAN MOLLMAN

Remarks:

Form 106-I8000