

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-003060 Instrument Location Walter J	City PD
Reason for Install/Repair	
Install After Receiving From Crime Laboratory [] Install After Location Change	
Other (Specify)	
Check When Done:	10.00
1. Surge Protector Installed/Property Grounded.	
2. Telephone Line Connected to Intoxilyzer® 8000.	
☑ 3. Breath Tube Heated.	
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio D	uring Test).
7. Repair and/or Maintenance Performed (if any):	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the	Intoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Ag	ency.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair (104-G), Print Test, ACA Test, and RFL Test.	Checkout (SFN59281, Form
Field Inspector Signature	Date 9/26/25
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determ the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Cricounty of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	
Reviewed/Certified By Anna Narehood	Certified Date 01Oct2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003060
Location = WCPD 8164.14.00 09/16
09/26/2025 15:56

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: Display Prelim Rslt? DABACABA Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? Standard Cyl #? 19825080A2 27 Standard Expiration? 09/05/2027 Oper No? 131581

Flow Cal. Date: 04/08/2008 Slope 702 Intercept -647610

IR Calibration Date: 06/04/2008

3um 9um

Oth Coef(*100): -5190 -7737

1st Coef(*100): 270580 138779

 0th Coef(*100):
 -5190
 -773

 1st Coef(*100):
 270580
 1387

 2nd Coef(*100):
 2499
 956

 H2O adj(mg/l*10k):
 377
 188

******* Printer Test End ********

Operator Signature
JOSHUA MCGOVEN

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WCPD
09/26/2025

Alcohol Analyzer SN 80-003060 8164.14.00 09/16 16:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:09
02 Std. Gas	0.079	16:09
03 Room Air	0.000	16:10
04 Std. Gas	0.079	16:10
05 Room Air	0.000	16:11
06 Std. Gas	0.079	16:11
07 Room Air	0.000	16:12

Lot No = 19825080A2

Cyl No = 27

Exp Date = 09/05/2027

County = 27

Oper No. = 131581

Operator Signature
JOSHUA MCGOVEN

Remarks: (1 Good test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WCPD
09/26/2025

SN 80-003060 8164.14.00 09/16 16:12

DRY CAL CHECK

Test			AC	Time
01	Room	Air	RFI*	16:13
02	Room	Air	RFI*	16:13

*Invalid Test Inhibited - RFI

Lot No = 19825080A2

Cyl No = 27

Exp Date = 09/05/2027

County = 27

Oper No. = 131581

1160

Operator Signature JOSHUA MCGOVEN

Remarks:

Form 106-I8000